

Hippocampal Moderate Sparing Enhanced Coverage Whole Brain (HMS-EC-WB) Model Description

Purpose:

This document describes the context in which the Hippocampal Moderate Sparing Enhanced Coverage Whole Brain (HMS-EC-WB) Model should be used, as well as how it was configured and validated. All instruction needed to use this model in your clinic can be found when you **read the first five pages of this document**.

Applicability:

- Note** *RapidPlan knowledge-based planning and its models are not intended to replace clinical decisions, provide medical advice or endorse any particular radiation plan or treatment procedure. The patients' medical professionals are solely responsible for and must rely on their professional clinical judgment when deciding how to plan and provide radiation therapy.*
- Note** *The performance of the HMS-EC-WB model may vary depending on the contouring and planning guidelines. Each site should validate the model with institution-specific contouring and planning guidelines before clinical use.*
- Note** *You should validate every DVH estimation model before using it clinically. This applies to any model, whether Varian provided, peer provided or the models you create yourself.*

- This model is designed to be used for RapidArc treatment plans for whole brain utilizing a hippocampal sparing technique. This model was created using 4 arc Halcyon plans calculated with AcurosXB but has also been validated with multiple (coplanar and non-coplanar) beam geometries on TrueBeam. (See Annex A for quantification of performance for each method).
- **This HMS-EC-WB model differs from prior HSWBv2 model in that this model reduces sparing to the Hippocampus to increase PTV Rx coverage from 95% to ~98.5%. Unless you are sure the intent is to trade hippocampal sparing for better Rx coverage near the hippocampus structure, it is recommended considering, instead, the previously released HSWBv2 RapidPlan model (or try both).**

Intent	<u>HLS-EC-WB</u> Limited Sparing Enhanced Coverage 20Gy/5fx-only	<u>HMS-EC-WB</u> Moderate Sparing Enhanced Coverage 30Gy Scalable Rx	<u>HSWBv2</u> Aggressive Sparing 30Gy Scalable Rx
Hippocampus Dmin	7.6Gy (38% Rx) ALARA	9Gy (30% Rx) ALARA	9Gy (30% Rx) ALARA
Hippocampus D0.03cc	13Gy (65% Rx)	16Gy (53.3% Rx)	16Gy (53.3% Rx) ALARA
PTV Rx dose coverage	20Gy @ 98%-99%	30Gy @ 98%-99%	30Gy @ ≥95%

- The model is intended to be used in conjunction with a **MU objective** with a **strength of 80** and **minimum MU 1000 (1600 for 6X-FFF)** and **maximum MU of 2500**. This MU objective **must be added manually each time**.
- The **“Automatic Intermediate Dose”** function of the Photon Optimizer was utilized with **MR3** return and **convergence mode: extended** selected in the calculation options. These **settings should be changed prior to starting the optimization or plan quality will be compromised**. Also, to reach desired homogeneity goals, consider an **additional intermediate dose optimization: “2xMR3”** (See Annex A)

- The model is intended for whole brain with hippocampal sparing without a simultaneous integrated boost (SIB) to gross disease. If SIB is intended to be utilized for boost volume, clinical validation of model performance is necessary. The model was not generated or fully validated for SIB clinical cases (See Annex D).

Target and OAR contouring and planning guidelines:

The HMS-EC-WB model was created using the following guidelines. Every patient must have a planning CT. The CT simulation scan must encompass the entire head to include the most superior aspect of the patient through the entire head. Axial slice thickness should not exceed 2.5mm and smaller axial cuts are recommended. The use of MRI guided contouring is also recommended. The MRI axial slice thickness should match the CT slice thickness as much as possible. It is recommended to obtain gadolinium-enhanced studies to include three-dimensional spoiled gradient (SPGR), magnetization-prepared rapid gradient echo (MP-RAGE), or turbo field echo (TRF) axial MRI scans with axial slice thickness not greater than 1.5mm. Standard axial and coronal gadolinium contrast-enhanced T1-weighted sequence and axial T2/FLAIR sequence scans should be acquired with no greater than 2.5mm slice thickness.

The planning target volumes (PTV) and the organs at risk (OARs) are contoured on the planning CT.

Target contouring guidelines:


Target name	Guidelines
CTV (or Brain)	Whole brain parenchyma through the foramen magnum
PTV_WB	CTV + 2mm setup margin added in this model, but 0-3mm are valid (institutional preference) excluding the hippocampal avoidance region

OAR contouring guidelines:

OAR name	Guidelines
Hippocampus(R+L)	Bilateral hippocampal contours (contoured in one structure); will be generated from the CT simulation image set fused to MRI image sets <u>Average Total Volume of this structure was 4.5cc</u> in training set cases. The largest volumes seen in the training set were ~7cc (>7.5cc were excluded).
Brainstem	Best generated on MRI image set and verified on CT simulation image set. Inferior aspect at the level of foramen magnum (should be at the interface of CTV inferior aspect) and superiorly to include midbrain
Spinal Canal	Superior aspect to begin at the distal edge of the brainstem through inferior aspect of the image set
Lens(R/L)	Use CT image set only for creation; bilateral contoured separately
Optic Nerve (R/L)	Use CT image set only for creation; bilateral contoured separately
Eye(R/L)	Best generated from CT simulation image set; delineate the entire globe of the eye; bilateral contoured separately
Optic Chiasm	Structure best visualized on MRI image set and confirmed on CT image set; located above the pituitary fossa which is located within the sella turcica
Lacrimal Gland(R/L)	Structure sits superior and lateral to the globe of the eye; best generated on CT simulation image set and verified on MRI; bilateral contoured separately

Optimization structure Guidelines*

Hippocampus+05(R+L)	Avoidance used to create PTV_WB for evaluation only (not trained in this model)
PTV_WBopt04	PTV excluding hippocampi +4mm additional margin from hippocampus structure. Reduced margin in optimization to achieve Rx @ 99% PTV_WB).
_Brain&BODY (Ring)	+20mm from the brain, SUB +5mm from brain, removed from outside body
_Brainstem#Hi (BS_PTV)	Brainstem SUB Hippocampus + 5mm
_Eyes&BODY (Face)	Eyes +100mm, SUB brain +20mm, removed from outside of the body



* Optimization structures can be automatically created with ESAPI PlanScoreCard tool (Annex C3)

Treatment planning guidelines:

All cases used to train and to validate the model were planned using head-first supine position with head positioned in a neutral position. All patients were immobilized with an aquaplast mask. A four arc VMAT technique was utilized with four full coplanar arcs on Halcyon with MLC in SX2 mode (standard mode on all Halcyon D / Drive and above configurations). Arcs had alternating clock-wise and counter clock-wise gantry rotations with collimator positions set at 315, 0, 45 and 90. The coplanar arcs had 359.8 degrees of arc rotation for each field. Arcs were positioned at a single isocenter located in the center of the target.

Full validation with different number of arcs, geometries and dose calculation methods on TrueBeam (M120 MLC) can be seen in Annex A, including a quantification of the relative dosimetric performance of each method.

The following dose prescription and planning guidelines were used for the cases to train and validate the model.

Target	PTV_WB	30Gy in 10 fractions (scales to alternate Rx, including 20Gy in 5 fractions)
	coverage	D100% at 99%; D98% > 30Gy; D2% < 31.5Gy (normalization to D100% >=99%)

OARs	Chiasm	D0.03cc < 31.5Gy
	Brainstem	D0.03cc < 33Gy
	Cord	D0.03cc < 33Gy
	Optic Nerve	D0.03cc < 31.5Gy
	Eye	Mean dose < 3Gy; Max dose <16.5Gy
	Lacrimal Gland	Mean dose < 4Gy
	Lens	D0.03cc < 4Gy
	Hippocampus	D0.03cc <16Gy; Mean dose < 9Gy; D100% < 9Gy (As Low As Reasonably Achievable)

References for contouring and planning guidelines:

Roberge D, Chan M, Gondi V. **CCTG CE. 7: Stereotactic Radiosurgery Compared With Hippocampal-Avoidant Whole Brain Radiotherapy (HA-WBRT) Plus Memantine for 5 or More Brain Metastases**

<https://www.ctg.queensu.ca/public/brain/brain-disease-site>

(HMS-EC-WB-2023 <https://medicalaffairs.varian.com/wholebrain-moderate-hippocampalsparing-30gy-vmat2>)

Liu H, Clark R, Magliari A, Foster R, Reynoso F, Schmidt M, Gondi V, Abraham C, Curry H, Kupelian P, Khuntia D, Beriwal S. **RapidPlan hippocampal sparing whole brain model version 2-how far can we reduce the dose?** Med Dosim. 2022 Autumn;47(3):258-263. doi: 10.1016/j.meddos.2022.04.003.

[https://www.meddos.org/article/S0958-3947\(22\)00039-5/fulltext](https://www.meddos.org/article/S0958-3947(22)00039-5/fulltext)

(HSWBv2-2022 <https://medicalaffairs.varian.com/wholebrain-hippocampalsparing-vmat2>)

Magliari V, Magliari A, Foster R. **Hippocampal Sparing Whole Brain: Rapid Plan Model Following the NRG-CC001 Protocol.** AAMD Conf Poster Present.

<https://medicalaffairs.varian.com/download/PosterPresentationAAMD2017RapidplanHCSWB.pdf>

(HSWBv1-2016 <http://medicalaffairs.varian.com/wholebrain-hippocampalsparing-vmat1>)

Brown P, Gondi V **NRG-CC001: A Randomized Phase III Trial of Memantine and Whole-Brain Radiotherapy With or Without Hippocampal Avoidance in Patients with Brain Metastases** <https://www.nrgoncology.org/Clinical-Trials/NRG-CC001>

Prokic V, et al **Whole Brain Irradiation with Hippocampal Sparing and Dose Escalation on Multiple Brain Metastases: A Planning Study on Treatment Concepts** <http://dx.doi.org/10.1016/j.ijrobp.2012.02.036>

Physicians with considerable experience in treating patients under the **CCTG CE.7 protocol**, additional planning goals and contouring consistency were completed to their clinical preference.

Structure codes:


To ensure robust structure matching between new cases and the structures defined in the model, it is recommended to use the following structure code assignment:

Structure name example	Structure name in model	Structure code(s) in model
PTV_WBopt04mm	PTV_WBopt	(PTV_High,PTV_Int Target)
Hippocampus_Totl	Hippocampus(R+L)	(275020)
Chiasm	OpticChiasm	(62045)
Brainstem	Brainstem	(79876)
Spinal Cord	SpinalCanal	(9680, 7647)
LOptic	OpticNerve(R/L)	(50878, 50875)
ROptic	OpticNerve(R/L)	(50878, 50875)
LEye	Eye(R/L)	(12515, 125124)
REye	Eye(R/L)	(12515, 125124)
LLacrimonal	Lacrimonal(L/R)	(59103, 59102)
RLacrimonal	Lacrimonal(L/R)	(59103, 59102)
LLens	Lens((R/L)	(58243, 58242)
RLens	Lens((R/L)	(58243, 58242)
NS_Ring	_Brain&BODY	(Control Region)
BrainstemPTV	_Brainstem#Hi	(Control Region)
Face	_Eyes&Body	(Control Region)

Optimization objectives and settings:

The following optimization objectives were defined in the model and will be generated when the model is applied to a new case:

Applying the model will also set the following parameters for the NTO:



MU objective is also recommended for planning with the following parameters:

Minimum 1600 for 6X-FFF, 1000 for 6X

******* MUST be added manually *******

Target	ID	Vol [%]	Dose	Priority	gEUD a
Yes	_BRAINSTEM#HI (Control)				
	Lower	99.9	100.5 %		210
Yes	PTV_WBopt4 'TV_High, PTV_Intermediate)				
	Upper	0.0	105.6 %		400
	Upper	1.5	103.5 %		210
	Upper	12.0	102.4 %		195
	Lower	99.5	101.6 %		280
	Lower	99.7	100.5 %		300
	Lower	99.9	95.7 %		220
	Upper gEUD		101.8 %		135 40.0
	_BRAIN&BODY (Control)				
	Upper	0.0	99.0 %		350
	_Eyes&BODY (Control)				
	Mean		Generated		85
	Line (preferring OAR)		Generated	Generated	110
	Brainstem (79876)				
	Upper	0.0	101.0 %		200
	Eye(R/L) (12515, 12514)				
	Upper (fixed vol., generated dose)	0.0	Generated		125
	Mean		Generated		95
	Line (preferring target)		Generated	Generated	85
	Hippocampus(R+L) (275020)				
	Upper	0.0	15,000 Gy		350
	Upper (fixed vol., generated dose)	99.0	Generated		190
	Upper (fixed vol., generated dose)	0.0	Generated		125
	Mean		Generated		100
	Line (preferring target)		Generated	Generated	95
	Lacrimonal(L/R) (59103, 59102)				
	Mean		Generated		65
	Line (preferring target)		Generated	Generated	80
	Lens(R/L) (58243, 58242)				
	Upper (fixed vol., generated dose)	0.0	Generated		110
	Mean		Generated		68
	Line (preferring target)		Generated	Generated	74
	OpticChiasm (62045)				
	Upper	0.0	100.5 %		160
	OpticNerve(R/L) (50878, 50875)				
	Upper	0.0	100.5 %		160
	Line (preferring target)		Generated	Generated	30
	SpinalCanal (9680, 7647)				
	Upper	0.0	100.5 %		160

Model Training:

This Hippocampal Moderate Sparing Enhanced Coverage Whole Brain (HMS-EC-WB) model was trained with the same final 42 case multi-institution CT dataset from the HSWBv2, structures were modified as needed. Each case was simulated with aquaplast mask immobilization and neutral head position.

All cases were initially re-planned to 30Gy in 3Gy fractions with 6X-FFF energy on a Varian Halcyon with SX2 MLC mode.

All cases utilized VMAT technique. Arcs had alternating clockwise and counterclockwise gantry rotations with collimator positions set at 315, 0, 45 and 90 degrees. The coplanar arcs had 359.8 degrees of arc rotation and were positioned with isocenter located in the center of the target.

The recursive method of model creation was utilized to generate a RapidPlan model with very consistent, high-quality plans developed with tight DVH prediction bands allowing for finely balanced hippocampal sparing, target coverage and homogeneity optimization objectives to be used. HMS-EC-WB uses plans created from HSWBv2 (which was, in-turn, created from HSWBv1 model released in 2016) as its starting point. These initial starting point plans were created without the HSWBv2 hippocampal DVH prediction line objectives and instead DVH point objectives were generated along the hippocampus line objectives. These objectives were then offset by fixed percentages toward higher dose levels to account for the dose gradient shifting toward the hippocampal structures in order to achieve the desired target coverage goal (Rx dose covering 99% of PTV_WB). A modified version of the V2.0 scorecard was created which adjusted previous metrics and added additional metrics to capture CE. 7 and it's author's clinical preference (aggressive target coverage and maximal homogeneity). The free PlanScoreCard ESAPI scoring tool was not only utilized to score plans (scores which guide tuning the model's automatically generated optimization objective priorities), but was also used to automatically create optimization structures (see Annex C3).

Those initial plans created from HSWBv2 model with offset hippocampal sparing objectives became the training set for the initial HMS-EC-WB model. A recursive model creation process was employed to ensure the final HMS-EC-WB training set consisted, exclusively, of plans generated from the initial HMS-EC-WB model. Evaluating plan scores at each step in the process informed multiple iterations of re-tuning the optimization objective set.

However, it was discovered through extensive objective tuning attempts on the final model training set that the offset hippocampal objectives from the plans for the initial training set cases were not quite aggressive enough. This resulted in prediction bands and auto created optimization objectives that were not adequate to spare the hippocampus to the degree desired on the scorecard. Therefore, hippocampal sparing objectives were purposely overshoot during optimization on the cases used to train the final model (even to the detriment of their score) to achieve higher scoring plans from with the auto created optimization objectives (with no manual offsets needed) using the final model. See table on the next page for details.

Model Validation:

The HMS-EC-WB model was validated using the 42 cases included in the final model training set. See the table on the next page to understand the model creation / validation and scores achieved throughout the process.

Five additional validation cases not included in the model training set are explored, in detail, in the Annex. Those cases are scored with various arc geometries (both coplanar and non-coplanar) on both Halcyon and TrueBeam (Millennium120 MLC) with differing numbers of intermediate dose optimizations, with a quantification of the relative dosimetric performance (score) of each method. Further validation including differing beam energies, dose calculation methods and convergence mode options, are available in the clinical description included with the HSWBv2 Rapidplan model.

V2.0 Scorecard=142		HMS-EC-WB Scorecard = 158.5 Total points		
Final Result		created by modV2.0 model with offset objectives	created by initial HMS-EC- WB with further offset obj*	Final Result
Patient	Final V2.0 Model	Training Set for initial HMS-EC-WB	Training Set for final HMS-EC-WB*	Final HMS-EC- WB Model
Patient 1	132.35	137.17		140.14
Patient 3	132.47	136.37		134.52
Patient 4	137.17	146		146.26
Patient 5	129.81	140.69		141.48
Patient 8	132.53	136.73		136.49
Patient 9	132.6	135.36		135.6
Patient 11	132.32	131.73		133.95
Patient 13	131.92	135.29		140.56
Patient 16	136.6	146.44		146.22
Patient 19	134.64	139.24		138.92
Patient 20	131.66	132.07		132.13
Patient 21	131.36	133.81		134.26
Patient 23	134.14	139.89		140.72
Patient 24	133.82	141.03		141.03
Patient 25	134.88	139.53		137.89
Patient 27	133.96	141.26		141.62
Patient 28	131.32	133.53		131.95
Patient 30	133.23	138.34		140.85
Patient 34	131.16	135.78		136.87
Patient 35	132.72	135.71		136.14
Patient 44	129.91	134.27		125.48
Patient 45	129.64	132.47		134.31
Patient 47	133.63	141.32		142.91
Patient 48	135.3	143.18		143.28
Patient 49	134.7	143.56		144.12
Patient 50	131.7	137.23		138.78
Patient 52	132.01	134.59		136.65
Patient 54	134.41	140.28		141.77
Patient 55	132.49	133.25		133.47
Patient 57	134.69	139.7		140.43
Patient 60	131.74	133.29		133.25
Patient 64	134.42	138.93		139.33
Patient 65	133.45	136		139.27
Patient 66	132.58	130.75		131.95
Patient 68	132.24	138.88		138.64
Patient 69	130.8	124.41		126.85
Patient 70	133.5	134.59		134.95
Patient 71	132.54	135.08		134.27
Patient 72	132.72	135.85		135.08
Patient 77	131.67	135.78		137.46
Patient 80	134.76	141.45		143.19
Patient 85	135.24	140.6		140.3
Average	132.9714286	137.18	*	137.70

*plans generated from initial model used further manual offset hippocampal objectives for sufficiently aggressive DVH prediction and auto-created optimization objectives to enhance final model's performance, lower scores

Annex Directory

Annex A: Validation Results

A1 Beam Arrangements: Halcyon and TrueBeam

A2 1xMR3, 2xMR3, 3xMR3 (Convergence Mode: Extended)

A3 Rapidplan model v15.6 and v18 builds and optimization + calculation versions

Annex B: Comparison to alternate Hippocampal Sparing RapidPlan models: coverage/sparing

B1 30Gy Isodose DVH & Score comparison HMS-EC-WB/HSWBv2 @ 30Gy

B2 Scorecard differences between alternate RapidPlan models

B3 20Gy Isodose DVH & Score comparison HLS-EC-WB/HMS-EC-WB/HSWBv2 @ 20Gy

Annex C: Scorecard Details

C1 HMS-EC-WB full Scorecard

C2 PlanScoreCard ESAPI tool: where to find

C3 PlanScoreCard ESAPI tool: automatically generate derived structures

Annex D: Examples applying this model for Simultaneous Integrated Boost

D1 Example SIB Plan 27.5Gy in 5Fx

D2 SIB DVH Comparison

D3 Planning Structures

D4 Metastasis proximity to hippocampus

Annex E: Acknowledgements

Annex F: Distribution and compatibility

Annex A: Validation Results

A1: Beam Arrangements (6X-FFF, AcurosXB v17, extended convergence mode, MR3 return, 2x Intermediate dose)

HMS-EC-WB additional 5 case validation on Halcyon and TrueBeam (M120 MLC)

Halcyon SX2MLC: coplanar collimator: 315°, 0°, 45°, 90°.

Truebeam M120MLC

HyperArc: full 4 arc arrangement

4 Arcs Non-Coplanar: 2 full arcs 0° couch 315°/45°



collimators and two vertex 180° (PA) -> 5° (from AP)

90° couch CW/CCW paired arcs with 315°/45° collimator

3 Arcs Coplanar: collimator: 315°, 45°, 90°


4 Arcs Coplanar: collimator same as 3 Arcs except

90° split X jaw superior/inferior to hippocampus



Patient	Halcyon	TrueBeam			
	4 Arcs (Coplanar)	4 Arcs (Non-Coplanar)	3 Arcs (Coplanar)	4 Arcs (Coplanar)	HyperArc (Non-Coplanar)
36	131.85	134.09	125.96	124.73	132.19
37	140.76	134.81	131.33	132.92	131.66
39	137.85	135.40	129.64	130.98	134.99
40	136.94	133.87	126.48	126.55	133.04
41	134.07	134.00	126.96	126.81	131.35
Average	136.29	134.43	128.07	128.40	132.65

Patient 36 selected DVH:



Fields	Dose	Field Alignments	Clinical Goals	Optimization Objectives	Dose Statistics	Reference Points	Calculation Models	Plan Sum	Dose Cover [%]	Sampling Cover [%]	Min Dose [Gy]	Max Dose [Gy]	Mean Dose [Gy]
Hippocampi	Unapproved	3FTB	HMS_EC_RP18CM17	3.2	100.0	100.0	8.722	16.427	12.235				
PTV_3000	Unapproved	3FTB	HMS_EC_RP18CM17	1228.7	100.0	100.0	19.005	34.056	30.706				
Hippocampi	Unapproved	4FNC1B	HMS_EC_RP18CM17	3.2	100.0	100.0	7.970	16.261	11.405				
PTV_3000	Unapproved	4FNC1B	HMS_EC_RP18CM17	1228.7	100.0	100.0	20.392	33.806	30.794				
Hippocampi	Unapproved	4FTB	HMS_EC_RP18CM17	3.2	100.0	100.0	8.839	16.404	12.237				
PTV_3000	Unapproved	4FTB	HMS_EC_RP18CM17	1228.7	100.0	100.0	17.447	34.343	30.719				
Hippocampi	Unapproved	HALL4rcs1	HMS_EC_RP18CM17	3.2	100.0	100.0	7.754	16.017	11.631				
PTV_3000	Unapproved	HALL4rcs1	HMS_EC_RP18CM17	1228.7	100.0	100.0	20.345	33.794	30.747				


Annex A: Validation Results

A2: 1xMR3, 2xMR3, 3xMR3 (Convergence Mode: Extended)


(v17 algorithms, extended convergence mode, MR3 return)

HMS-EC-WB patient 36 validation on Halcyon

Always use extended convergence mode and MR3 return:




1XMR3= "Automatic Intermediate dose"




*GPU not required

2XMR3 and 3XMR3: multiple Intermediate dose optimizations, "current plan as an intermediate dose for optimization"



Patient 36 Scores / DVH:

·HCSWB-036	HAL4Arcs1xMR3	129.69	2249 MU
·HCSWB-036	HAL4Arcs2xMR3	131.85	2560 MU
·HCSWB-036	HAL4Arcs3xMR3	133.20	2626 MU



Annex A: Validation Results


A3: Rapidplan v15.6 and v18 model versions (Halcyon 4 arc, 6X-FFF, AcurosXB v17 and v18 with ELM)

This model was created on a V18 Eclipse system with v17 Photon Optimizer and AcurosXB. For backwards compatibility, all training set cases were exported and imported into a V15.6 system. The V15.6 model was trained and the optimization objectives from the V18 model were replicated. Finally, the V15.6 model was exported from the V15.6 system back into the V18 system and plans were reoptimized with scores comparing with the results of the V18 native model with both PO & AcurosXB v17 and PO & AcurosXB v18 with Enhanced Leaf Modeling.

	Halcyon		
	4 Arcs (Coplanar)		
Patient	RPv18 PO+AXBv17	RPv15.6 PO+AXBv17	RPv15.6 PO+AXBv18wELM
36	131.85	133.8	127.78
37	140.76	136.82	142.21
39	137.85	136.52	139.58
40	136.94	133.75	137.09
41	134.07	135.05	130.15
Average	136.29	135.12	135.36

Annex B: Comparison of HMS-EC-WB / HSWBv2: different tradeoffs (coverage vs sparing)

B1 Isodose, DVH, score comparison: HMS-EC-WB & HSWBv2models at 30Gy (patient 36)



Plan Score Comparison 30Gy

RP model:	V2.0 Scorecard (142 points)				HMS-EC-WB Scorecard (158.5 points)			
	HSWBv2		HMS-EC-WB		HSWBv2		HMS-EC-WB	
patient 36	132.08	93.01%	111.06	78.21%	125.18	78.98%	131.85	83.19%
patient 37	133.24	93.83%	112.64	79.32%	127.46	80.42%	140.76	88.81%
patient 39	132.17	93.08%	115.12	81.07%	123.15	77.70%	137.85	86.97%
patient 40	133.39	93.94%	115.53	81.36%	123.35	77.82%	136.94	86.40%
patient 41	131.82	92.83%	118.49	83.44%	123.42	77.87%	134.07	84.59%
average	132.54	93.34%	114.568	80.68%	124.512	78.56%	136.294	85.99%


The HSWBv2 RapidPlan model, released in 2022, and its associated dosimetric scorecard express a precise clinical intent: maximal hippocampal sparing with WB-PTV coverage at $\geq 95\%$, among other quality metrics. In 2023, alternative hippocampal sparing RapidPlan models were released along with a specific dosimetric scorecard which quantifies a plan's ability reach a new desired intent. Here we compare two plans, both meeting the NRG-CC01 protocol's basic dosimetric goals, but beyond those goals further specific dosimetric intent can be realized. In this case, either by further sparing of the hippocampi (HSWBv2) or via enhanced prescription coverage by constraining the low dose into only the hippocampal structures and enforcing full Rx dose coverage by the end of the hippocampal avoidance zone (hippocampus+05). This 2023 RapidPlan model variation is thereby named "Hippocampal Moderate Sparing Enhanced Coverage Whole Brain" or HMS-EC-WB. The relevant dosimetric scorecard for each RapidPlan model scores the dose made by its associated model and the non-associated model to demonstrate the utility of the dosimetric scorecards to capture a precise clinical intent and quantify the relative success in implementing a dose plan to realize said intent. More of the total points are achieved by the single click RapidPlan model which aims for the specific dosimetric intent laid out in each model's associated scorecard as demonstrated on the five model validation cases in the table (dose and DVH shown only for patient 36).

B2 Clinical intent with precision: HLS-EC-WB, HMS-EC-WB, HSWBv2 Scorecards Overview


Id	StructureId	Metric Text	Max Score	Id	StructureId	Metric Text	Max Score	Id	StructureId	Metric Text	Max Score
0	PTV_2000	Volume at 20Gy [%]	20	0	PTV_3000	Volume at 30Gy [%]	20	0	PTV_3000	Volume at 30Gy [%]	15
1	PTV_2000	Dose at 98% [Gy]	14	1	PTV_3000	Dose at 98% [Gy]	14	1	PTV_3000	Dose at 98% [Gy]	14
2	PTV_2000	Dose at 95% [Gy]	2	2	PTV_3000	Dose at 95% [Gy]	2	2	PTV_3000	Dose at 2% [Gy]	11
3	PTV_2000	Dose at 2% [Gy]	11	3	PTV_3000	Dose at 2% [Gy]	11	3	PTV_3000	Volume at 105% [%]	5.5
4	PTV_2000	Volume at 105% [%]	5.5	4	PTV_3000	Volume at 105% [%]	5	4	PTV_3000	Dose at 0.03CC [Gy]	4.5
5	PTV_2000	Dose at 0.03CC [Gy]	8	5	PTV_3000	Dose at 0.03CC [Gy]	2.5	5	PTV_3000	HI [1 - 99]/30]	2
6	PTV_2000	MaxDose [Gy]	5	6	PTV_3000	MaxDose [Gy]	2.5	6	PTV_3000	Conformation No. at [28.5Gy]	1
7	PTV_2000	HI [1 - 99]/20]	2	7	PTV_3000	HI [1 - 99]/30]	2	7	Hippocampi	Dose at 0.03CC [Gy]	7.5
8	PTV_2000	Conformation No. at [12.7Gy]	1	8	PTV_3000	Conformation No. at [28.5Gy]	1	8	Hippocampi	MeanDose [Gy]	12
9	Hippocampi	Dose at 0.03CC [Gy]	7.5	9	Hippocampi	Dose at 0.03CC [Gy]	12	9	Hippocampi	Dose at 100% [Gy]	17
10	Hippocampi	MeanDose [Gy]	12	10	Hippocampi	MeanDose [Gy]	16	10	OpticChiasm	Dose at 0.03CC [Gy]	3.5
11	Hippocampi	Dose at 100% [Gy]	17	11	Hippocampi	Dose at 100% [Gy]	17	11	BrainStem	Dose at 0.03CC [Gy]	3
12	OpticChiasm	Dose at 0.03CC [Gy]	3.5	12	OpticChiasm	Dose at 0.03CC [Gy]	3.5	12	SpinalCord	Dose at 0.03CC [Gy]	3.5
13	BrainStem	Dose at 0.03CC [Gy]	3	13	BrainStem	Dose at 0.03CC [Gy]	3	13	OpticNerve_L	Dose at 0.03CC [Gy]	3.5
14	SpinalCord	Dose at 0.03CC [Gy]	3.5	14	SpinalCord	Dose at 0.03CC [Gy]	3.5	14	OpticNerve_R	Dose at 0.03CC [Gy]	3.5
15	OpticNerve_L	Dose at 0.03CC [Gy]	3.5	15	OpticNerve_L	Dose at 0.03CC [Gy]	3.5	15	Eye_L	MaxDose [Gy]	2
16	OpticNerve_R	Dose at 0.03CC [Gy]	3.5	16	OpticNerve_R	Dose at 0.03CC [Gy]	3.5	16	Eye_L	MeanDose [Gy]	2
17	Eye_L	MaxDose [Gy]	2	17	Eye_L	MaxDose [Gy]	2	17	Eye_R	MaxDose [Gy]	2
18	Eye_L	MeanDose [Gy]	2	18	Eye_L	MeanDose [Gy]	2	18	Eye_R	MeanDose [Gy]	2
19	Eye_R	MaxDose [Gy]	2	19	Eye_R	MaxDose [Gy]	2	19	LacrimalGland_L	MeanDose [Gy]	3.5
20	Eye_R	MeanDose [Gy]	2	20	Eye_R	MeanDose [Gy]	2	20	LacrimalGland_R	MeanDose [Gy]	3.5
21	LacrimalGland_L	MeanDose [Gy]	3.5	21	LacrimalGland_L	MeanDose [Gy]	3.5	21	Lens_L	Dose at 0.03CC [Gy]	2.25
22	LacrimalGland_R	MeanDose [Gy]	3.5	22	LacrimalGland_R	MeanDose [Gy]	3.5	22	Lens_R	Dose at 0.03CC [Gy]	2.25
23	Lens_L	Dose at 0.03CC [Gy]	2.25	23	Lens_L	Dose at 0.03CC [Gy]	2.25	23	_Brain&BODY	Volume at 99.5% [CC]	5
24	Lens_R	Dose at 0.03CC [Gy]	2.25	24	Lens_R	Dose at 0.03CC [Gy]	2.25	24	_Brain&BODY	MaxDose [%]	5
25	_Brain&BODY	Volume at 99.5% [CC]	5	25	_Brain&BODY	Volume at 99.5% [CC]	5	25	_BrainStem#Hi	Dose at 95% [Gy]	2
26	_Brain&BODY	MaxDose [%]	5	26	_Brain&BODY	MaxDose [%]	5	26	_Eyes&BODY	MeanDose [Gy]	5
27	_BrainStem#Hi	Dose at 95% [Gy]	2	27	_BrainStem#Hi	Dose at 95% [Gy]	2	27			
28	_Eyes&BODY	MeanDose [Gy]	5	28	_Eyes&BODY	MeanDose [Gy]	5	28			
HLS-EC-WB Scorecard Total			158.5	HMS-EC-WB Scorecard Total			158.5	HSWBv2 Scorecard Total			142

The above table is an overview of the differences between the dosimetric scorecard point distribution for the various clinical intents each variation represents. Hippocampal Limited Sparing (HLS-EC-WB) on the left has more points allocated to reducing hotspots in the target whereas Moderate Sparing (HMS-EC-WB) allocates those points to further hippocampal sparing. These can both be contrasted with the older 2022 scorecard with fewer metrics and total points. The above table does not include the specific piecewise linear scoring functions and the specific scoring value ranges where each metric fails (0 points) and where max points are awarded. The point distribution differences between HMS-EC-WB and its neighboring scorecards are highlighted. See Annex C1 for full example HMS-EC-WB dosimetric scorecard.

B3 Isodose, DVH, score: HLS-EC-WB, HMS-EC-WB(scaled to 20Gy), HSWBv2(scaled to 20Gy)



Plan Score Comparison 20Gy (V2.0 and Moderate scorecard and plan doses scaled *.667)																		
model	HLS-EC-WB Scorecard (158.5 points)				HMS-EC-WB Scorecard (158.5 points)				V2.0 Scorecard (142 points)									
	HSWBv2		HMS-EC-WB		HLS-EC-WB		HSWBv2		HMS-EC-WB		HLS-EC-WB		HSWBv2		HMS-EC-WB		HLS-EC-WB	
patient 36	123.42	77.87%	135.43	85.44%	141.32	89.16%	125.18	78.98%	131.85	83.19%	124.67	78.66%	132.08	93.01%	111.06	78.21%	110.48	77.80%
patient 37	128.49	81.07%	139.63	88.09%	143.93	90.81%	127.46	80.42%	140.76	88.81%	119.22	75.22%	133.24	93.83%	112.64	79.32%	112.35	79.12%
patient 39	116.72	73.64%	134.61	84.93%	143.31	90.42%	123.15	77.70%	137.85	86.97%	133.58	84.28%	132.17	93.08%	115.12	81.07%	117.58	82.80%
patient 40	114.66	72.34%	129.7	81.83%	137.3	86.62%	123.35	77.82%	136.94	86.40%	120.13	75.79%	133.39	93.94%	115.53	81.36%	116.68	82.17%
patient 41	123.76	78.08%	130.27	82.19%	141.54	89.30%	123.42	77.87%	134.07	84.59%	121.98	76.96%	131.82	92.83%	118.49	83.44%	112.31	79.09%
average	121.41	76.60%	133.928	84.50%	141.48	89.26%	124.512	78.56%	136.294	85.99%	123.916	78.18%	132.54	93.34%	114.568	80.68%	113.88	80.20%



Dose	Reference Points	Dose Statistics
Structure	Approval Status	Plan
Hippocampi	Unapproved	HLS-EC-WB
PTV_2000	Unapproved	HLS-EC-WB
Hippocampi	Unapproved	HMS-EC-WBx20G
PTV_3000	Unapproved	HMS-EC-WBx20G
Hippocampi	Unapproved	V2.0-20Gy
PTV_3000	Unapproved	V2.0-20Gy

Structure	Volume [cm ³]	Dose Cover [%]	Sampling Cover [%]	Min Dose [Gy]	Max Dose [Gy]	Mean Dose [Gy]
Hippocampi	3.2	100.0	100.2	5.022	13.287	8.215
PTV_2000	1369.3	100.0	100.0	15.594	22.376	20.589
Hippocampi	3.2	100.0	100.2	5.215	10.773	7.823
PTV_3000	1228.7	100.0	100.0	13.684	22.689	20.680
Hippocampi	3.2	100.0	100.2	3.884	8.755	5.325
PTV_3000	1228.7	100.0	100.0	8.014	22.703	20.447

The HSWBv2 RapidPlan model, released in 2022, is configured to be scalable in Rx dose and works well to create 20Gy plans too. In 2023, alternative hippocampal sparing RapidPlan models were released along with their relevant specific dosimetric scorecards which quantify a plan's ability reach the desired intent. Here we compare three plans, all meet more relaxed sparing goals of an institution specific 20Gy protocol, but beyond those goals further specific dosimetric intent can be realized. In 2023, an additional RapidPlan model variation was released named "Hippocampal Limited Sparing Enhanced Coverage Whole Brain". HLS-EC-WB 20Gy does not well control the max dose in the hippocampus (and will likely fail that NRC-CC01 dose goal if scaled to 30Gy) but instead maximizes dose homogeneity. The relevant dosimetric scorecard for each RapidPlan model (scaled to the same 20Gy Rx) scores the dose made by its associated model and the non-associated models to demonstrate the utility of the dosimetric scorecards to capture a precise clinical intent and quantify the relative success in implementing a dose plan to realize said intent. More of the total points are achieved by the single click RapidPlan model designed to implement the specific dosimetric intent laid out in each model's associated scorecard as demonstrated on the five model validation cases in the table (dose and DVH shown only for patient 36).

Annex C: Dosimetric scorecard details and PlansScoreCard ESAPI tool

C1 Full Scorecard for HMS-EC-WB: Validation patient 36



Download full Scorecard(json), DICOM case example and this RapidPlan model:

<https://medicalaffairs.varian.com/wholebrain-moderate-hippocampalsparing-30gy-vmat2>

Annex C: Dosimetric scorecard and PlansScoreCard ESAPI tool

C2 PlanScoreCard ESAPI tool: where to find

[Varian-MedicalAffairsAppliedSolutions \(https://github.com/Varian-MedicalAffairsAppliedSolutions/MAAS-PlanScoreCard\)](https://github.com/Varian-MedicalAffairsAppliedSolutions/MAAS-PlanScoreCard)

Varian-MedicalAffairsAppliedSolutions / MAAS-PlanScoreCard Public

Notifications Fork 8 Star 13

<> Code Issues 5 Pull requests Discussions Actions Projects Security Insights

main 5 branches 43 tags Go to file Code

About

Medical Affairs Applied Solutions ESAPI tool to create ScoreCards and score plans; in-metric Boolean/expansion; normalize dose to max score; multi-patient batch scoring

Readme View license 13 stars 10 watching 8 forks Report repository

Releases 4

V16.1-PlanScoreCard-V3.1.7.12-0... (Latest) on Apr 3 + 3 releases

Packages

No packages published

Contributors 6

Languages

C# 100.0%

varian-ma Update README with batch mode screen and typo 345bc08 on May 4 253 commits

.github/workflows	Update V18 action to use latest v18 ESAPI package	5 months ago
NormalizeToScorecard	Testing to run normalization through the same application, but those ...	2 years ago
PlanScoreCard	Added commented option for resolving expiration	2 months ago
.gitattributes	Add .gitignore and .gitattributes.	2 years ago
.gitignore	Add .gitignore and .gitattributes.	2 years ago
BasicInstallQuickStart.md	Update BasicInstallQuickStart.md	10 months ago
ChangeLog.md	Update ChangeLog.md	3 months ago
FAQ.md	Update FAQ.md	6 months ago
InstallGuidePart2IntoSystemScriptsDi...	Rename InstallGuidePart2IntoSystemScriptsDirectory.md to InstallGuide...	10 months ago
PlanScoreCard.sln	Update github actions to fix missing release attachment problem and e...	5 months ago
README.md	Update README with batch mode screen and typo	last month
Troubleshooting.md	Create Troubleshooting.md	10 months ago
license.txt	added license.txt	last year

☰ README.md

PlanScoreCard

Medical Affairs Applied Solutions ESAPI tool to create dosimetric ScoreCards and score plans.

Features:

- Quantitative piecewise linear scoring functions for each metric
 - optional: flag for point where "variation acceptable" sited on referenced protocol
 - optional: note section to site referenced protocol or justification for metric (points)
 - optional: qualitative colors and labels for metric points, ie: orange="Just OK"
- Advanced scoring criteria supported
 - ConformationNumber
 - ConformityIndex
 - DoseAtVolume
 - HomogeneityIndex

Currently, the source code is shared on the Varian Innovation Center GitHub where it can be downloaded and compiled with Visual Studio 2022 (including with the free community edition), now in the releases section users can find precompiled binaries ready to run in all compatible versions of Eclipse (v15.6+). PlanScoreCard is made available under the Varian Limited Use Software License Agreement.

Annex C: Dosimetric scorecard and PlansScoreCard ESAPI tool

C3 PlanScoreCard ESAPI tool: automatically generate derived structures

The PlanScoreCard tool has a feature where derived structures (made with Boolean and expansion tools) can be created automatically. These structures can be created temporarily (to be used for scoring the plan but never saved back to the database) or (if the ESAPI tool has been approved for writing) the PlanScoreCard tool's configuration file can be edited so these generated structures are saved.

Below are screen captures showing how to build structures

Ring Structure Generation

+20mm from the brain, SUB +5mm from the brain, and removed from outside of the body

Operation	Structure	Margin [mm]	+
SUB	Brain	- 20	+
SUB	Brain	- 5	+
AND	BODY	- 0	+


Group1
<Brain>[20 SUB <Brain>]5 AND <BODY>

Structure Builder Comment
<Brain>[20 SUB <Brain>]5 AND <BODY>

Structure Id:

FINALIZE STRUCTURE

Ring Structure




Annex C: Dosimetric scorecard and PlansScoreCard ESAPI tool


C3 PlanScoreCard ESAPI tool: automatically generate derived structures

Brainstem Target Structure Generation

Brainstem SUB Hippocampi + 5mm



Brainstem Target Structure



Annex C: Dosimetric scorecard and PlansScoreCard ESAPI tool

C3 PlanScoreCard ESAPI tool: automatically generate derived structures

Face Structure Generation

Eyes + 100mm, SUB Brain +20mm, and removed from outside of the body

Operation	Structure	Margin (mm)
AND	Leye	100
OR	Reye	100
SUB	Brain	20
AND	BODY	0


Group1
<Leye>[100 OR <Reye>[100 SUB <Brain>[20 AND <BODY>]

Structure Builder Comment
<Leye>[100 OR <Reye>[100 SUB <Brain>[20 AND <BODY>]

Structure Id:

FINALIZE STRUCTURE

Face Structure



Annex D: Examples applying this model for Simultaneous Integrated Boost

The following is taken from the related Limited Sparing 20Gy model (HLS-EC-20Gy) but is applicable for this model, especially when scaling this Moderate Sparing (HMS-EC-30Gy scalable) model to 20Gy. However, when boosting metastasis, it may be better to use the 2022 HSWBv2 model and scaling it to the desired whole brain dose, as these "Enhanced Coverage" models are meant to not "miss" metastasis near the hippocampus. When boosting metastasis, consider taking advantage of the more aggressive hippocampal sparing in the 2022 HSWBv2 model, instead.


This model was trained for only the PTV whole brain target. However, it could be used to create SIB plans by cropping the PTV_WBopt, with some additional margin, from the PTVBoost target(s). The PTV whole brain should also be removed from the high risk PTV + 7mm to evaluate heterogeneity within the target. In the below examples, an additional 7mm margin was also used between the WB_PTVopt (20Gy) and the PTV_Boost (27.5Gy) target.

After cropping additional margins out of the PTV whole brain and PTV_WBopt, use HMS-EC-WB to automatically populate the optimizer as intended. Manually add upper and lower dose constraints for the PTVBoost, per prescription. Ensure that the margin removed from the WB_PTVopt, to accommodate the PTVBoost, is not too conservative or aggressive. Consider increasing the MU objective Maximum MU to >2500 if utilizing this RapidPlan model for SIB treatment plans.

When using this method, the HMS-EC-WB model has no knowledge of the higher dose level target when generating the DVH prediction bands and relative optimization objectives. This situation could cause the model to create objectives that are no longer relevant for your patient and could, as a result, create undesirable plans. The degree to which the objectives are off relates directly to the distance that the higher dose target is from the OARs. This is especially important for the hippocampus DVH bands being predicted and the increased dose the high risk PTV is prescribed relative to the PTV whole brain. In scenarios where the PTVboost is near the hippocampi, it is advised to copy the hippocampi into an evaluation structure and that is cropped with an additional margin away from the PTVboost. This hippocampi evaluation structure is to then be matched to the hippocampi in the HMS-EC-WB model for DVH estimation and optimization. Due to these various clinical scenarios, the usage of SIB cannot be endorsed by the creators of this model. However, what each user does with this model is at the discretion of the user and their associated clinical, physics, and medical staff.


D1: Example SIB Plan 27.5Gy in 5Fx

PTV_2750_Total boost to 27.5Gy and PTV_2000_SIB to 20Gy (Halcyon 4 arc, 6X-FFF, AcurosXB v17).



D2: SIB DVH Comparison


DVH comparison: SIB (HCSWB_SIB) vs Non-SIB (HCSWB HMS-EC-WB) plans



D3: SIB planning structures


PTV_2000 structure cropped 5mm from hippocampus and 7mm from PTV_2750

PTV_OPT optimization structure cropped 4mm from hippocampus and 7mm from PTV_2750




D4: PTV boost proximity to hippocampi

CCTG CE. 7 allows for hippocampal sparing wherever possible, even with metastasis close to or overlapping. In the below example, to achieve desired hippocampal sparing with this proximity to the boost volume, the Hippocampus is copied and cropped by 5mm from the PTV boost. This hippocampus_Eval structure is then matched to the hippocampus structure in the HMS-EC-WB model for optimization and evaluation.



D5: Metastases proximity to hippocampi (sequential boost)

The following workflow can be used in the occurrence for when contoured brain metastases are in proximity of the hippocampus and sequential boosting is implemented. To maintain prescription dose (20Gy) to the contoured GTV and PTV brain metastases, the hippocampus is copied and cropped by 5mm from the "PTV met". This hippocampus_Eval structure is then later matched to the hippocampus structure in the HMS-EC-WB model for optimization and evaluation. After the PTV_2000 structure and PTV_OPT structures are cropped 5mm and 4mm from the hippocampus respectively, the PTV Met with an additional 2mm margin can be added back to the both the PTV_2000 and PTV_OPT. This allows the model to account for the desired gradient to achieve coverage of the PTV met and reduce dose to the hippocampus_Eval.



Annex E: Acknowledgements

Manually optimized plans created by Ryan Clark, MS CMD

All data generated and compiled by Ryan Clark, MS CMD, Lesley Rosa, CMD and Anthony Magliari, MS CMD

Dosimetric Scorecard HMS-EC-WB and derived structures designed by Anthony Magliari, MS CMD and Ryan Clark, MS CMD

Model generated plans created by Lesley Rosa, CMD, Ryan Clark, MS CMD and Anthony Magliari, MS CMD

Clinical Description document created by Anthony Magliari, MS CMD, Ryan Clark, MS CMD

Please reference future publication:

Dosimetric Scorecards express clinical intent with precision: alternate hippocampal sparing whole brain RapidPlan models with enhanced target coverage, homogeneity and less aggressive hippocampal sparing

Kareem Rayn, Anthony Magliari, Ryan Clark, Lesley Rosa, Robert Doucet, Line Comeau, Alan Nichol, Russel Ruo, David Roberge

Annex F: Distribution and compatibility

This RapidPlan model is to be distributed exclusively via the links found on Varian Medical Affairs:

<https://medicalaffairs.varian.com/wholebrain-moderate-hippocampalsparing-30gy-vmat2>

Please do not re-distribute this model as number of downloads will be tracked (strictly to judge the success of this project).

This RapidPlan model was created, tested, and rebuilt with both Eclipse v18.0 and v15.6. For older versions of Eclipse (v13.x), please find the older HSWBv1.