



Proton FLASH Ultra High Dose Rate Treatment Plans: Dosimetric Comparisons to Standard of Care

Anthony Magliari, MS, CMD
2023 AAMD National Meeting

FLASH therapy is under development and not available for commercial sale.

I am employed by Varian on the
Medical Affairs team

My job currently includes testing new
products and providing feedback
including FLASH radiotherapy

I've been creating FLASH treatment
plans since 2017

I used development builds of Eclipse
and Non Clinical modes of the
ProBeam delivery system for
examples

FLASH therapy is under development and not available for commercial sale.

DISCLAIMER:

varian

A Siemens Healthineers Company

"The views expressed in this presentation are mine, and mine alone. They do not represent those of Varian, A Siemens Healthineers Company"

Previous AAMD national meeting talks:
2016 Role of the Dosimetrist in the Knowledge Based
Planning Era
2017 IMRT and VMAT: current and future best practices
2019 FLASH Radiotherapy: A Look at Ultra-high Dose Rate
Research and Treatment Plans
2022 Knowledge Sharing and the Power of Plan Challenge
ScoreCards


Disclaimer 2: FLASH Intellectual Property

- Varian intends to leverage all intellectual property protections and is committed to protecting our innovations throughout the world.
- Varian has numerous pending and issued patents worldwide related to FLASH therapy and technology

FLASH therapy is under development and not available for commercial sale.

- 1 FLASH Review
- 2 Defining Dose Rate
- 3 Dosimetric Comparisons
- 4 Future Possibilities

LET'S
DO
THIS



FLASH therapy is under development and not available for commercial sale.




FLASH Review

Background slides to explain Flash Radiotherapy


1

FLASH therapy is under development and not available for commercial sale.

Ultra High Dose Rates



Electron Flash reduces lung fibrosis in mice compared to conventional electron therapy



**NORMAL
TISSUE**

**SIGNIFICANT
FIBROSIS**

**NO APPARENT
DAMAGE**

**Control
0 Gy/s**


**CONV 17 Gy
0.03 Gy/s**

**Flash 17 Gy
60 Gy/s**

Favaudon et al., Science Trans Med 2014


FLASH therapy is under development and not available for commercial sale.

The promise of Flash




FLASH therapy is under development and not available for commercial sale.


The promise of Flash



The promise of Flash




Consistent evidence of tumor control and reduction in toxicity—across anatomical sites



Thorax, Leg, and Lung research conducted with mice. Brain research conducted with rats.

FLASH therapy is under development and not available for commercial sale.

The FLASH effect is a biological effect




Significance

- UHDR reduces toxicities in animal models
- UHDR remains efficient to control tumor in animal models
- The combination of these two outcomes is named the "FLASH effect"

Vozenin et al, Clinical cancer research, 2019

FLASH therapy is under development and not available for commercial sale.

The first exhaustive dose escalation study comparing CONV and FLASH in mice



Significance

In mice, the calculated pFLASH sparing factor is around 40%.

Concretely and considering single dose regimen, the FLASH factor allowed a 10-15Gy dose escalation

Sorensen et al., Radiother & Oncol, 2022

FLASH therapy is under development and not available for commercial sale.

FLASH Therapy – a hot topic in cancer treatment

>200 peer-reviewed articles

Varian and the Cincinnati Children's/UC Health Proton Therapy Center Announce Initial Patient Treated in the FAST-01 First Human Clinical Trial of FLASH Therapy for Cancer

Oncology
November 19, 2020

PALO ALTO, Calif., and CINCINNATI, Ohio, Nov. 19, 2020 /PRNewswire/ -- Varian (NYSE: VAR) and the Cincinnati Children's/UC Health Proton Therapy Center today announce the start of the first clinical trial of FLASH therapy as part of the recently opened FAST-01 study (FeAsibility Study of FLASH Radiotherapy for the Treatment of Symptomatic Bone Metastases). The clinical trial involves the investigational use of Varian's ProBeam® particle accelerator modified to enable radiation therapy delivery at ultra-high dose rates (dose delivered in less than 1 second) and is being conducted at the Cincinnati Children's/UC Health Proton Therapy Center with John C. Breneman M.D., Medical Director of the center, serving as principal investigator.

The FLASH-enabled ProBeam® system is an investigational device and is limited by United States law for investigational use.



Defining Dose Rate

Pencil Beam Scanning Dose Rate


2

FLASH therapy is under development and not available for commercial sale.

External Beam Radiation

Plane with point of interest


Radiation Source



FLASH therapy is under development and not available for commercial sale.

Beam's Eye View


Plane with point of interest



FLASH therapy is under development and not available for commercial sale.

Conventional Delivery

Proton
(Spot Scanning)




Photon/Electron
(Pulsed Broad Beam)




FLASH therapy is under development and not available for commercial sale.

Conventional Dose Delivery

Proton
(Spot Scanning)




Photon/Electron
(Pulsed Broad Beam)




FLASH therapy is under development and not available for commercial sale.

FLASH Dose Delivery (1/10th speed)

Proton
(Spot Scanning)




Photon/Electron
(Pulsed Broad Beam)




FLASH therapy is under development and not available for commercial sale.

Average Dose Rate

Photon/Electron FLASH




- Dose Rate = $\frac{\Delta \text{Dose}}{\Delta \text{Time}} = 50 \text{ Gy/s}$
- “field” dose rate




FLASH therapy is under development and not available for commercial sale.

Average Dose Rate

Proton FLASH




- Field Dose Rate = $\frac{\Delta \text{Dose}}{\Delta \text{Time}} = 50 \text{ Gy/s}$
- Conventional




FLASH therapy is under development and not available for commercial sale.

Pencil Beam Scanning Dose Rate

Proton FLASH



- “PBS” Dose Rate = $\frac{\Delta \text{Dose}}{\Delta \text{Time}} = 130 \text{ Gy/s}$



FLASH therapy is under development and not available for commercial sale.

Dose Rate Peer-Reviewed Article

Published September 2020




<https://doi.org/10.1002/mp.14456>

FLASH therapy is under development and not available for commercial sale.

Dose rate and dose accumulation in PBS field


Dose accumulation over time depends on the location within the treatment field and the spot scanning pattern



FLASH therapy is under development and not available for commercial sale.

Dose rate and dose accumulation in PBS field


Dose accumulation over time depends on the location within the treatment field and the spot scanning pattern



FLASH therapy is under development and not available for commercial sale.

Dose and Dose Rate Display



Single Field in a multi-field Flash Lung Treatment Plan



FLASH therapy is under development and not available for commercial sale.

Dose and Dose Rate Line Profile


Through Normal Lung Tissue



FLASH therapy is under development and not available for commercial sale.

Rethinking Treatment Planning

How to Review a Flash Plan to ensure the Flash conditions are met?




FLASH therapy is under development and not available for commercial sale.

Transmission planning

How to deliver at UHDR?

- Need
Ultra-High Dose Rate should be available for treatment to get in Flash range.
- Challenge
Dose rate drops significantly due to beam current loss.
- Solution
No use of upstream energy degrader.
Highest dose rate achieved at Maximum Beam Energy from cyclotron: **250 MeV**
First approach: Bragg Peak is outside of Patient



FLASH therapy is under development and not available for commercial sale.

Optimization for Flash

How to optimize on both dose and dose rate?

- Need

Increase dose rate without compromising the dose distribution

- Challenge

Dose distribution suffers when dose rate (MUs) increase due


- Solution:

- Optimize using the Min MU as a surrogate for dose rate (high min MU = high current = high dose rate)



Increasing min MU to increase dose rate

Drastic degradation of plan dosimetry...if not done correctly





Dose Rate Optimization

Standard Dose Volume Objectives:

minimum MU objective + priority

Structure	Constraint	Type	Priority	Volume (%)	Value	
Esophagus	DVC	max	50	0.59	3.00	Gy
Esophagus	DVC	max	50	5.22	1.00	Gy
Esophagus	DVC	max	50	2.83	2.09	Gy
Esophagus	DVC	max	50	7.58	0.25	Gy
Esophagus	DVC	max	75	0.00	3.50	Gy
SpinalCord	DVC	max	75	0.00	3.00	Gy
SpinalCord	DVC	max	75	1.49	2.00	Gy
SpinalCord	DVC	max	75	2.89	1.00	Gy
PTV	DVC	max	110	0.00	35.00	Gy
PTV	DVC	min	250	100.00	34.00	Gy
Lungs-GTV	DVC	max	50	1.59	8.37	Gy
Lungs-GTV	DVC	max	50	6.81	3.33	Gy

minimum MU: 400 Priority: 400




FLASH therapy is under development and not available for commercial sale.

Optimization for Flash

How to optimize on both dose and dose rate?

- Need
Maintain a homogenous dose distribution while increasing dose rate
- Challenge
Fewer spots with higher MUs per spot can lead to uneven dose distribution (hot or cold spots)
- Solution:
Place spots evenly to ensure homogenous dose distribution without hot spots



FLASH therapy is under development and not available for commercial sale.

Concept Review:

Best case scenario for achieving ultra-high dose rate Flash plans when utilizing the PBS dose rate definition:

- a. Large targets with a low dose per fraction
- b. Large targets with a high dose per fraction
- c. Small targets with a low dose per fraction
- d. Small targets with a high dose per fraction



Concept Review:

Proton spot spacing and layer spacing's effect on treatment quality is:

- a. Inconsequential, using the defaults as defined in the treatment planning system is best
- b. It matters and is best to have as many spots as possible with dense spot and layer spacing for the best dosimetric plan quality
- c. Spending time finding an optimal sparse spacing of spots can lead to better dosimetry and larger minimum MU per field, also decreasing the treatment time






Dosimetric Comparison

FLASH plans compared to IMPT standard of care


3

FLASH therapy is under development and not available for commercial sale.

Understanding SOBP vs 250MeV broad-beam vs 250MeV single spot




Understanding SOBP vs 250MeV broad-beam vs 250MeV single spot




Prostate SBRT with Urethral Sparing: 5 field 250MeV vs 2 field IMPT

PTV37.5Gy SIB Prostate 40Gy, 5fx 10Gy/fx



Prostate SBRT with Urethral Sparing: 5 field 250MeV vs 2 field IMPT

PTV37.5Gy SIB Prostate 40Gy, 5fx 10Gy/fx



Structure	Reference Points	Dose Statistics	Plan	Score	Value	Score	Max	Metric Plot
Urethra	Unapproved	V200POW0	5field	2.3	100.0	0.425	100.0	28.41%
Urethra	Unapproved	SChex-bl-init	2field-IMPT	37.6	100.0	100.0	3.893	35.47%
Rectum	Unapproved	V200POW0	5field	37.6	100.0	0.000	0.000	49.12%
Rectum	Unapproved	SChex-bl-init	2field-IMPT	37.0	100.0	0.000	0.000	40.96%
PTV	Unapproved	V200POW0	5field	84.2	100.0	38.279	48.000	47.38%
PTV	Unapproved	SChex-bl-init	2field-IMPT	84.2	100.0	24.401	48.000	42.29%
Femurs	Unapproved	V200POW0	5field	100.0	100.0	0.159	0.159	9.65%
Femurs	Unapproved	SChex-bl-init	2field-IMPT	100.0	100.0	1.420	23.115	16.53%
Bladder	Unapproved	V200POW0	5field	100.4	100.0	0.000	0.000	43.08%
Bladder	Unapproved	SChex-bl-init	2field-IMPT	100.4	100.0	0.000	0.000	44.08%
Urethra	Unapproved	V200POW0	5field	1.0	100.0	0.000	0.000	0.00%

FLASH therapy is under development and not available for commercial sale.

Prostate SBRT with Urethral Sparing: 5 field 250MeV vs 2 field IMPT

PTV37.5Gy SIB Prostate 40Gy, 5fx 10Gy/fx

[N-Opt] V200POW0100.5: 150.74/220.00 (68.52%)
[2field-IMPT] SChex-bl-init: 104.57/220.00 (47.53%)


Structure	Metric Text	Max Score
PTV	Volume at 40Gy [%]	5
PTV	Volume at 36.25Gy [%]	5
PTV	Dose at 0.07% [Gy]	8
PTV-CTV	Dose at 0.15% [Gy]	3
PTV-UrethraPRV	Volume at 36.25Gy [%]	15
PTV-UrethraPRV	Volume at 40Gy [%]	20
PTV-Urethra	Dose at 99.99% [Gy]	10
PTV-Urethra	Dose at 99.99% [Gy]	5
PTV-Urethra	Volume at 30Gy [%]	10
CTV-UrethraPRV	Volume at 40Gy [%]	10
CTV-UrethraPRV	Dose at 99.99% [Gy]	12
Rectum	Volume at 36.25Gy [%]	4
Rectum	Volume at 42Gy [%]	3
Rectum	Volume at 34.43Gy [%]	3
Rectum	Dose at 50% [Gy]	8
Rectum	Dose at 20% [Gy]	6
Rectum	Dose at 0.5% [Gy]	5
Bowel_Small	Volume at 35Gy [%]	6
Bowel_Small	Volume at 25Gy [%]	4
Bowel_Large	Volume at 35Gy [%]	3
Bowel_Large	Volume at 25Gy [%]	2
Bladder	Volume at 36.25Gy [%]	4
Bladder	Volume at 42Gy [%]	4
Bladder	Volume at 32.62Gy [%]	6
Bladder	Dose at 40% [Gy]	5
Urethra	Volume at 35Gy [%]	10
Urethra	Dose at 0.1% [Gy]	10
Femurs	Volume at 19.5Gy [%]	3
Fiducials	Dose at 1% [Gy]	1
PTV_RING	Dose at 0.01% [Gy]	5
BODY_RING	Dose at 0.01% [Gy]	5

Plan Scores:	RapidPlan: [N-Opt] V200POW0100.5: 150.74/220.00 (68.52%)	RapidPlan: [2field-IMPT] SChex-bl-init: 104.57/220.00 (47.53%)					
Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV	Volume at 40Gy [%]	V200POW0100.5	85.44 %	0.35	5.00	
			SChex-bl-init	82.85 %	0.00	5.00	
		GOOD[4] IDEAL[5]					
2	PTV	Volume at 36.25Gy [%]	V200POW0100.5	94.32 %	2.43	5.00	
			SChex-bl-init	92.81 %	1.58	5.00	
		GOOD[4.5] IDEAL[5]					
3	PTV	Dose at 0.07% [Gy]	V200POW0100.5	44.62 Gy	4.49	8.00	
			SChex-bl-init	46.92 Gy	0.00	8.00	
		IDEAL[8] GOOD[6.5] VARIATION[0]					
4	PTV-CTV	Dose at 0.15% [Gy]	V200POW0100.5	43.97 Gy	0.04	3.00	
			SChex-bl-init	45.25 Gy	0.00	3.00	
		IDEAL[3] GOOD[2.5] VARIATION[0]					
5	PTV-UrethraPRV	Volume at 36.25Gy [%]	V200POW0100.5	99.86 %	14.86	15.00	
			SChex-bl-init	98.82 %	13.82	15.00	
		VARIATION[0] GOOD[11] IDEAL[15]					

FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx

5 field transmission flash



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx[®]

ROG 0915

Table 1: Conformity of Prescribed Dose for Calculations Based on Deposition of Photon Beam Energy in Heterogeneous Tissue


PTV Volume (cc)	Ratio of Prescription Isodose Volume to the PTV Volume		Ratio of 50% Prescription Isodose Volume to the PTV Volume, R _{50%}		Maximum Dose (in % of dose prescribed) @ 2 cm from PTV in Any Direction, D _{2cm} (Gy)		Percent of Lung Receiving 20 Gy Total or More, V ₂₀ (%)	
	Deviation		Deviation		Deviation		Deviation	
	None	Minor	None	Minor	None	Minor	None	Minor
1.8	<1.2	<1.5	<5.9	<7.5	<50.0	<57.0	<10	<15
3.8	<1.2	<1.5	<5.5	<6.5	<50.0	<57.0	<10	<15
7.4	<1.2	<1.5	<5.1	<6.0	<50.0	<58.0	<10	<15
13.2	<1.2	<1.5	<4.7	<5.8	<50.0	<58.0	<10	<15
22.0	<1.2	<1.5	<4.5	<5.5	<54.0	<63.0	<10	<15
34.0	<1.2	<1.5	<4.3	<5.3	<58.0	<68.0	<10	<15
50.0	<1.2	<1.5	<4.0	<5.0	<62.0	<77.0	<10	<15
70.0	<1.2	<1.5	<3.5	<4.8	<66.0	<86.0	<10	<15
95.0	<1.2	<1.5	<3.3	<4.4	<70.0	<89.0	<10	<15
126.0	<1.2	<1.5	<3.1	<4.0	<73.0	>91.0	<10	<15
163.0	<1.2	<1.5	<2.9	<3.7	<77.0	>94.0	<10	<15

Note 1: For values of PTV dimension or volume not specified, linear interpolation between table entries is required.

Note 2: Protocol deviations greater than listed here as "minor" will be classified as "major" for protocol compliance (see Section 6.7).

FLASH therapy is under development and not available for commercial sale.


SBRT Lung 34Gy 1fx ScoreCard Example (143.5 total points)



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx


5 field transmission flash | CI: 1.33 | Plan Score: 136.46



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx


5 field transmission flash MinMU=600 40Gy/sec Dose Rate Threshold DVH



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx


2 field IMPT | CI: 1.39 | Plan Score: 136.03



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx


3 field IMPT | CI: 1.36 | Plan Score: 138.42



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx


12 field IMRT (HD-120) | CI: 1.08 | Plan Score: 136.56



FLASH therapy is under development and not available for commercial sale.

SBRT Lung 34Gy 1fx

Halcyon 4 partial arc | CI: 1.05 | Plan Score: 139.49



SBRT Lung 34Gy 1fx


2 field IMPT, 3 field IMPT, 5 field FLASH comparison

Lung PTV=23.69cc	Score x/143.5	CI
FLASH 5 fields	136.46	1.33
IMPT 2 fields	136.03	1.39
IMPT 3 fields	138.42	1.36
IMRT 12 fields	136.56	1.08
VMAT 4 partial arcs	139.49	1.05

FLASH therapy is under development and not available for commercial sale.

GBM 24Gy/36Gy SIB 6fx

2 field IMPT vs 5 field transmission flash



FLASH therapy is under development and not available for commercial sale.

GBM 24Gy/36Gy SIB 6fx

5 field transmission flash vs. 2 field IMPT

PLAN SCORECARD ***NOT VALIDATED FOR CLINICAL USE**

SCORE CARD SELECTION

Score Card ID: Plan Rx matches scorecard.
GBM11_36-24GySIB6fxJOK.pon

PLAN SELECTION

Patient ID: gbm11_FLASH

Patient	Course	Plan	Score (V23K)
gbm11_FLASH	IMPT	CSPTCSmR	188.54
gbm11_FLASH	IMPT	bik Final_6field	
gbm11_FLASH	IMPT	bik Plan1	
gbm11_FLASH	IMPT	bik Plan2	
gbm11_FLASH	IMPT	bik Ran3-NUPC	
gbm11_FLASH	IMPT	bik 6field-test	
gbm11_FLASH	250MeV	SIB.5	
gbm11_FLASH	250MeV	4FidSIB.5	
gbm11_FLASH	250MeV	SCSIB.5	
gbm11_FLASH	250MeV	SIB.6	
gbm11_FLASH	250MeV	SIB.7	
gbm11_FLASH	250MeV	iC4FidSIB.7	
gbm11_FLASH	250MeV	4FidSIB.6	
gbm11_FLASH	250MeV	4FidSIB.7	
gbm11_FLASH	250MeV	SIB.8	227.19

Norm Export Print

Bound by the terms of the VarianUSLA

Plan Scores: gbm11_FLASH: [IMPT] SCPTCSmRS: 188.54/236.50 (79.72%)
gbm11_FLASH: [250MeV] SIB.8: 227.19/236.50 (96.06%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV24GyFlash	Dose at 99% [Gy]	SCPTCSmRS	22.20 Gy	99.00	99.00	[Plot]
			SIB.8	21.38 Gy	99.00	99.00	
2	PTV24GyFlashEval	Dose at 0.03CC [Gy]	SCPTCSmRS	37.43 Gy	0.09	5.00	[Plot]
			SIB.8	33.27 Gy	4.27	5.00	
3	PTV24GyFlashEval	Dose at 5CC [Gy]	SCPTCSmRS	34.32 Gy	2.28	5.00	[Plot]
			SIB.8	30.79 Gy	3.75	5.00	
4	PTV24GyFlash	Volume at 24Gy [%]	SCPTCSmRS	94.24 %	19.77	24.00	[Plot]
			SIB.8	96.27 %	24.00	24.00	
5	PTV36GyFlash	Volume at 36Gy [%]	SCPTCSmRS	40.48 %	0.00	24.00	[Plot]
			SIB.8	98.26 %	24.00	24.00	
6	PTV36GyFlash	Dose at 99% [Gy]	SCPTCSmRS	31.66 Gy	7.66	10.00	[Plot]
			SIB.8	35.72 Gy	10.00	10.00	
7	PTV36GyFlash	Dose at 0.03CC [Gy]	SCPTCSmRS	40.22 Gy	9.83	10.00	[Plot]
			SIB.8	40.19 Gy	0.87	10.00	

FLASH therapy is under development and not available for commercial sale.

Liver 50Gy 5fx (10Gy fractions) 5 field transmission flash vs VMAT

Metric Id	StructureId	Metric Text	Max Score
0	PTV Liver	Dose at 95% [Gy]	11
1	PTV Liver	Dose at 99% [Gy]	11
2	PTV Liver	Dose at 0.03CC [Gy]	11
3	Liver-GTV	MeanDose [Gy]	6
4	Liver-GTV	MeanDose [Gy]	5
5	Liver-GTV	Dose at 700CC [Gy]	10
6	Esophagus	Dose at 0.5CC [Gy]	6
7	Stomach	Dose at 0.5CC [Gy]	6
8	Duodenum	Dose at 0.5CC [Gy]	6
9	Bowel	Dose at 0.5CC [Gy]	6
10	Ring	Volume at 50Gy [%]	6
11	Ring	Dose at 0.05CC [Gy]	6
12	SpinalCordPRV	Dose at 0.5CC [Gy]	6
13	Kidneys	Dose at 33% [Gy]	6
14	Kidneys	MeanDose [Gy]	6
15	Kidneys	Dose at 10% [Gy]	6
16	Heart	Dose at 0.5CC [Gy]	6
17	Ribs	Dose at 2% [Gy]	6

Plan Scores: 2022RSS_Liver2: [5field] SCV300WO: 120.18/126.00 (95.38%)
2022RSS_Liver2: [VMAT] Hal1: 119.71/126.00 (95.01%)


Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV Liver	Dose at 95% [Gy]	SCV300WO	50.41 Gy	11.00	11.00	
			Hal1	49.11 Gy	10.64	10.64	
2	PTV Liver	Dose at 99% [Gy]	SCV300WO	50.03 Gy	11.00	11.00	
			Hal1	48.23 Gy	10.36	10.36	
3	PTV Liver	Dose at 0.03CC [Gy]	SCV300WO	52.61 Gy	10.99	11.00	
			Hal1	52.69 Gy	10.97	10.97	
4	Liver-GTV	MeanDose [Gy]	SCV300WO	8.59 Gy	5.85	6.00	
			Hal1	9.60 Gy	5.72	5.72	
5	Liver-GTV	MeanDose [Gy]	SCV300WO	8.59 Gy	3.35	5.00	
			Hal1	9.60 Gy	3.14	3.14	
6	Liver-GTV	Dose at 700CC [Gy]	SCV300WO	3.36 Gy	10.00	10.00	
			Hal1	6.96 Gy	9.35	9.35	
7	Esophagus	Dose at 0.5CC [Gy]	SCV300WO	9.93 Gy	5.67	6.00	
			Hal1	7.16 Gy	5.76	5.76	

[5field] SCV300WO: 120.18/126.00 (95.38%)
[VMAT] Hal1: 119.71/126.00 (95.01%)

FLASH therapy is under development and not available for commercial sale.

MultiMet test 3 mets on different planes, 1fx 21Gy

50%
10.5Gy




FLASH therapy is under development and not available for commercial sale.

MultiMet test

3 mets on different planes, 1fx 21Gy

10%
2.1Gy



FLASH therapy is under development and not available for commercial sale.

MultiMet test

3 mets on different planes, 1fx 21Gy

Metric Id	StructureId	Metric Text	Max Score
0	1PTV_RF_21	Volume at 21Gy [%]	40
1	2PTV_RC_21	Volume at 21Gy [%]	40
2	3PTV_LC_21	Volume at 21Gy [%]	40
3	1PTV_RF_21	Dose at 99.5% [Gy]	25
4	2PTV_RC_21	Dose at 99.5% [Gy]	25
5	3PTV_LC_21	Dose at 99.5% [Gy]	25
6	1PTV_RF_21	Dose at 0.03CC [Gy]	18
7	2PTV_RC_21	Dose at 0.03CC [Gy]	18
8	3PTV_LC_21	Dose at 0.03CC [Gy]	18
9	RingRx.1-1_PTV1	Dose at 0.03CC [Gy]	10
10	RingRx.1-1_PTV2	Dose at 0.03CC [Gy]	10
11	RingRx.1-1_PTV3	Dose at 0.03CC [Gy]	10
12	RingMid.3-1.5_PTV1	Volume at 10.5Gy [%]	5
13	RingMid.3-1.5_PTV2	Volume at 10.5Gy [%]	5
14	RingMid.3-1.5_PTV3	Volume at 10.5Gy [%]	5
15	RingMid.3-1.5_PTV1	Dose at 0.01CC [Gy]	20
16	RingMid.3-1.5_PTV2	Dose at 0.01CC [Gy]	20
17	RingMid.3-1.5_PTV3	Dose at 0.01CC [Gy]	20
18	RingLow1.0_PTVs	Volume at 6Gy [%]	5
19	RingLow1.0_PTVs	Volume at 4Gy [%]	5
20	RingLow1.0_PTVs	Volume at 2Gy [%]	5
21	Brainstem	Dose at 1CC [Gy]	5
22	Brainstem	Dose at 0.03CC [Gy]	2.5
23	Brainstem	Dose at 0.1CC [Gy]	2.5
24	Cochl_L	Dose at 0.03CC [Gy]	2.5
25	Cochl_L	Dose at 0.1CC [Gy]	2.5
26	Cochl_R	Dose at 0.03CC [Gy]	2.5
27	Cochl_R	Dose at 0.1CC [Gy]	2.5

Metric Id	StructureId	Metric Text	Max Score
28	Chiasm	Dose at 1% [Gy]	5
29	Chiasm	Dose at 0.1CC [Gy]	2.5
30	Chiasm	Dose at 0.03CC [Gy]	2.5
31	OptNrv_L	Dose at 1% [Gy]	5
32	OptNrv_L	Dose at 0.1CC [Gy]	2.5
33	OptNrv_L	Dose at 0.03CC [Gy]	2.5
34	OptNrv_R	Dose at 1% [Gy]	5
35	OptNrv_R	Dose at 0.1CC [Gy]	2.5
36	OptNrv_R	Dose at 0.03CC [Gy]	2.5
37	Eye_L	Dose at 0.03CC [Gy]	1
38	Eye_L	Dose at 0.1CC [Gy]	1
39	Eye_R	Dose at 0.03CC [Gy]	1
40	Eye_R	Dose at 0.1CC [Gy]	1
41	Lens_L	Dose at 0.03CC [Gy]	1
42	Lens_L	Dose at 0.1CC [Gy]	1
43	Lens_R	Dose at 0.03CC [Gy]	1
44	Lens_R	Dose at 0.1CC [Gy]	1


Plan Scores:	Ma_3metP4: [N-Opt] 5fieldWO: 337.31/427.00 (78.99%)	Ma_3metP4: [N-Opt] Halcyon: 377.51/427.00 (88.41%)					
Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	1PTVRF_21	Volume at 21Gy [%]	5fieldWO	99.50 %	40.00	40.00	
			Halcyon	99.50 %	40.00	40.00	
2	2PTVRC_21	Volume at 21Gy [%]	5fieldWO	99.84 %	40.00	40.00	
			Halcyon	99.73 %	40.00	40.00	
3	3PTVLC_21	Volume at 21Gy [%]	5fieldWO	99.67 %	40.00	40.00	
			Halcyon	99.91 %	40.00	40.00	
4	1PTVRF_21	Dose at 99.5% [Gy]	5fieldWO	20.99 Gy	24.75	25.00	
			Halcyon	21.00 Gy	25.00	25.00	
5	2PTVRC_21	Dose at 99.5% [Gy]	5fieldWO	21.32 Gy	25.00	25.00	
			Halcyon	21.00 Gy	25.00	25.00	

Ma_3metP4: [N-Opt] 5fieldWO: 337.31/427.00 (78.99%)
Ma_3metP4: [N-Opt] Halcyon: 377.51/427.00 (88.41%)

FLASH therapy is under development and not available for commercial sale.

SBRT spine 24Gy 1fx

10 field transmission flash vs. 3 field IMPT



FLASH therapy is under development and not available for commercial sale.

SBRT spine 24Gy 1fx

10 field transmission flash vs. 3 field IMPT

FLASH therapy is under development and not available for commercial sale.

Plan Scores:	[N-Opt] 10V75-99-600p: 277.31/310.00 (89.45%)	[N-Opt] SC3FLD-TEST9: 257.38/310.00 (83.02%)					
1	CTVEAL	Dose at 99.9% [Gy]	10V75-99-600p SC3FLD-TEST9	22.10 Gy 22.34 Gy	23.42 23.62	25.00 Score Stats Max=23.62 Mean=23.52 Min=23.42	Metric Plot
2	PTVEAL	Volume at 24Gy [%]	10V75-99-600p SC3FLD-TEST9	88.58 % 89.99 %	12.86 14.98	20.00 Score Stats Max=14.98 Mean=13.92 Min=12.86	Metric Plot
3	PTVEAL	Volume at 22.8Gy [%]	10V75-99-600p SC3FLD-TEST9	97.76 % 98.35 %	15.10 15.34	16.00 Score Stats Max=15.34 Mean=15.22 Min=15.10	Metric Plot
4	PTVEAL	Dose at 99.5% [Gy]	10V75-99-600p SC3FLD-TEST9	21.09 Gy 21.14 Gy	13.43 13.49	15.00 Score Stats Max=13.49 Mean=13.46 Min=13.43	Metric Plot
5	PTVEAL	Dose at 0.035CC [Gy]	10V75-99-600p SC3FLD-TEST9	28.80 Gy 30.25 Gy	13.00 5.15	15.00 Score Stats Max=13.49 Mean=9.07 Min=5.15	Metric Plot
6	PTVEAL	Volume at 28Gy [%]	10V75-99-600p SC3FLD-TEST9	1.53 % 2.05 %	11.86 11.12	14.00 Score Stats Max=11.86 Mean=11.49 Min=11.12	Metric Plot
7	PTVEAL	Volume at 26.2Gy [%]	10V75-99-600p SC3FLD-TEST9	27.32 % 22.20 %	7.26 8.37	10.00 Score Stats Max=8.37 Mean=7.82 Min=7.26	Metric Plot
8	CTV	Dose at 99.9% [Gy]	10V75-99-600p SC3FLD-TEST9	19.44 Gy 20.43 Gy	1.00 1.00	1.00 Score Stats Max=1.00 Mean=1.00 Min=1.00	Metric Plot
9	PTV	Volume at 24Gy [%]	10V75-99-600p	86.80 %	0.68	1.00	Metric Plot

SBRT spine 24Gy 1fx


10 field transmission flash vs. 3 field IMPT vs. failed options and non-FLASH

Number of fields	Min MU score assigned	Min MU target	Score /310
2 (ML-IMPT)	0	-	*
3 (ML-IMPT)	0	-	257.38
5 (250MeV)	0	-	*
10 (250MeV)	0	-	290.14
10 (250MeV)	40	150	287.82
10 (250MeV)	40	300	290.34
10 (250MeV)	40	450	290.28
10 (250MeV)	80	450	289.36
10 (250MeV)	r99	600	284.10
10 (250MeV)	r99*	600	277.31

FLASH therapy is under development and not available for commercial sale.

SBRT spine 24Gy 1fx

PBS Dose Rate analysis: Minimum MU=600 (Dose Rate Threshold DVH >40Gy/sec)



FLASH therapy is under development and not available for commercial sale.



Future Possibilities

Where FLASH proton planning could go next

4

FLASH therapy is under development and not available for commercial sale.

ScoreCard Based Optimization

Direct optimization on piecewise linear DVH based score function ranges

PLAN SCORECARD **NOT VALIDATED FOR CLINICAL USE**

SCORE CARD SELECTION

Score Card ID: SC_HeadAndNeck-AAMD2023_PlanStud

Plan Rx matches scorecard.

Plan Scores: AAMD2023PLANSTUDY: [C2] Vorono10MU: 64.87/150.00 (43.25%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV63	Volume at 63Gy [%]	Vorono10MU	100.00 %	10.00	10.00	
			Vorono10MU	99.93	9.90		
2	PTV60	Volume at 60Gy [%]	Vorono10MU	99.81 %	7.36	7.50	
			Vorono10MU	99.85	7.39		
3	PTV60	Volume at 63Gy [%]	Vorono10MU	98.49 %	0.00	2.50	
			Vorono10MU	0.62	2.47		
4	PTV57	Volume at 57Gy [%]	Vorono10MU	96.07 %	6.07	10.00	
			Vorono10MU	96.83	6.83		
5	PTV57	Volume at 59.85Gy [%]	Vorono10MU	92.76 %	0.00	2.50	
			Vorono10MU	30.92	0.95		
6	PTV54	Volume at 54Gy [%]	Vorono10MU	99.32 %	9.32	10.00	
			Vorono10MU	97.39	7.39		
7	PTV54	Volume at 56.7Gy [%]	Vorono10MU	98.31 %	0.00	2.50	

VarianPLAID5-02-2023 - PLAID for pFLASH

LAUNCH/QUIT IM-WEIGHT OPT POSITION OPT

BRAGG RM OPT SWAP SWAP-PP

INTERMEDIATE DOSE POST-PROCESS SCRIPT

optimization round complete
comment len (252 max): 194
Saving optimized plan ...
course id: C2
pland id: SCVorono10MU
NumberOfFractions: 30
prescribedDose: 63.0
planNo: 1.0
treatmentPercentage: 1.0
Applying optimized spot weights...
Calculating Dose...
Snapshot Saved!

Plot Controls
Total Score

Optimized Plan Score

RELOAD PLAN


Bound by the terms of the VarianUSLA

FLASH therapy is under development and not available for commercial sale.

ScoreCard Based Optimization

250MeV (non Flash) AAMD 2023 Plan Study Phase 2

AAMD2023PLANSTUDY: [C2] SCVoronoi10MU: 135.42/150.00 (90.28%)
AAMD2023PLANSTUDY: [C1] SCVoronoi10MU: 140.49/150.00 (93.66%)



Plan Name	Structure	Score Metric	Plan ID	Value	Score	Max	Min	Min/Max
1	PTV95	Volume at 50Gy (%)	SCVoronoi10MU	99.98 %	9.99	10.00		99.98
2	PTV90	Volume at 50Gy (%)	SCVoronoi10MU	99.80 %	9.83	10.00		99.80
3	PTV97	Volume at 50Gy (%)	SCVoronoi10MU	99.96 %	9.92	10.00		99.96
4	PTV94	Volume at 50Gy (%)	SCVoronoi10MU	47.97 %	4.67	47.97		47.97
5	CTV83	Dose at 95% (Gy)	SCVoronoi10MU	63.39 Gy	90.00	100.00		63.39
6	CTV80	Dose at 95% (Gy)	SCVoronoi10MU	60.27 Gy	90.00	100.00		60.27
7	CTV87	Dose at 95% (Gy)	SCVoronoi10MU	57.22 Gy	90.00	100.00		57.22
8	CTV84	Dose at 95% (Gy)	SCVoronoi10MU	54.33 Gy	90.00	100.00		54.33
9	CTV85	Volume of Region (0Gy) (%)	SCVoronoi10MU	0.00 %	5.00	5.00		0.00
10	SpinalCord	Dose at 0.0001 (Gy)	SCVoronoi10MU	38.75 Gy	4.38	5.00		38.75
11	Brainstem	Dose at 0.0001 (Gy)	SCVoronoi10MU	28.79 Gy	5.00	5.00		28.79
12	ParotidL	Volume at 30Gy (%)	SCVoronoi10MU	34.08 %	5.00	5.00		34.08
13	ParotidR	Volume at 30Gy (%)	SCVoronoi10MU	23.38 %	3.77	5.00		23.38
14	ParotidL	MeanDose (Gy)	SCVoronoi10MU	23.36 Gy	5.00	5.00		23.36


FLASH therapy is under development and not available for commercial sale.

FLASH therapy is under development and not available for commercial sale.

Spot Weight And Position Optimization

DIRECT MACHINE PARAMETER OPTIMIZATION FOR INTENSITY MODULATED PROTON THERAPY

P. Lansonneur*, A. Magliari, L. Rosa, J. Perez and M. Folkerts
 Varian, a Siemens Healthineers Company



Optimization of spot positions starting from a regular grid

Results




Fig. 2. Objective function evolution for conventional optimization and DMPO (case 1).

The evolution of the objective function during optimization is displayed in fig 2. The final cost reached with DMPO was two times lower than conventional optimization, with similar number of iterations. Unlike the conventional optimization scheme, all plans optimized with the DMPO algorithm passed the RTOG pre-specified metrics (see table 2).

Metric	Goal	Conventional	DMPO
PTV, V _{50Gy} (%)	> 90	99.0	99.0
PTV, D _{max} (Gy)	56.2	64.2	53.9
Ring80, D _{max} (Gy)	< 26.2	28.5	25.8
Ribs, D _{max} (Gy)	< 36.9	33.9	32.0
Ribs, D ₁₀₀ (Gy)	< 28.8	28.9	28.3
Spinal Cord, D _{max} (Gy)	< 18	8.0	7.6
Esophagus, D _{max} (Gy)	< 25.2	15.4	13.7
Lungs, D _{max} (Gy)	< 13.5	3.9	3.9

Table 2. Several dose metrics and associated passing criteria evaluated for conventional optimization and DMPO. Failing metrics are indicated in red, improved in orange and passing in green.

The dose distributions delivered above 40 Gy/ls are displayed for each case in figure 3. For all fields, most of the dose is delivered at a dose rate above 40 Gy/ls. Example field dose distributions are shown in figure 4.




Fig. 3. Dose delivered at dose rate above 40 Gy/ls.




Fig. 4. Field dose distributions as seen from the BEV for a beam optimized from a hexagonal grid and with DMPO.

Compared to conventional optimization on a regular grid, DMPO resulted in smaller "hot" regions of lower magnitude in the final dose distributions. The fraction of irradiated volume receiving at least 40 Gy/ls was above 91 % for lungs (excluding PTV) and above 92 % for esophagus (see table 3). Overall, DMPO resulted in significant plan quality improvement for all patients.

% FLASH	Esophagus	Lungs-PTV	Tracheobronchial
case 1	100 %	97 %	98 %
case 2	95 %	94 %	93 %
case 3	98 %	95 %	95 %
case 4	95 %	91 %	92 %

Table 3. Fraction of irradiated volume (dose > 2 Gy) above 40 Gy/ls several organs at risk.

Conclusion

We have proposed a new algorithm to optimize FLASH plans. Optimizing both the spot weights and positions leads to better plan quality than conventional (flux grid) optimization. This work will support the creation of fractionalized IMPT plans for FLASH-RT. The authors thank M. Rogo and P. Nemeis for their help. The DMPO algorithm is covered by a pending patent application.

References

- [1] G. Vidotto et al. RTOG-0915 (2012).
- [2] P. Lansonneur et al. PTCOG-60 poster (2022).
- [3] B. Nelms, Practical Radiation Oncology 2, 296-305 (2012).
- [4] M. Folkerts et al. Medical Physics (2020).


Bragg Peak Proton FLASH - AAPM 2023

FLASH therapy is under development and not available for commercial sale.


End-to-end verification of multi-field-optimization FLASH-IMPT using a clinical GBM case

Miriam Krieger¹, Pierre Lansonneur², Uli Weber³, Christoph Schuy³, Yuri Simeonov⁴, Petar Penchev⁴, Marta Rovituso⁵, Ernst van der Wal⁵, Anthony Magliari⁶, Christopher Bayer⁶, Aaron Marshall⁶, Gerard Paquet⁶, Klemens Zink⁴, Michael M Folkerts⁶

- 1 Varian Medical Systems Particle Therapy GmbH & Co
- 2 Varian Medical Systems, Le Plessis-Robinson, France
- 3 GSI Helmholtzzentrum für Schwerionenforschung, D
- 4 Technische Hochschule Mittelhessen, Giessen, Germ
- 5 HollandPTC, Delft, The Netherlands
- 6 Varian Medical Systems Inc., Palo Alto, CA, USA



FLASH therapy is under development and not available for commercial sale.



Thank you!

Special Thanks to:

- Michael Folkerts
- Jessica Perez
- Pierre Lansonneur
- Miriam Krieger
- Lesley Rosa



FLASH therapy is under development and not available for commercial sale.