

Taking Knowledge Based Planning to the Next Level - Modern Tools to Build Better Models Faster

Ryan Clark, MS, CMD. R.T.(T)
AAMD Annual Meeting 2023

1

Disclosure

- ▶ I am employed by Varian
- ▶ The views expressed in this presentation are mine and mine alone and do not represent those of Varian.

2

2

Outline

- ▶ Review of basic Knowledge Based Planning concepts
 - Challenges of the H&N treatment site
- ▶ Background for this H&N model
 - AI clustering by parotid involvement
 - Specific model for “Cluster 0” (both parotids partially involved)
 - Case selection for training set cases within cluster 0 (3 targets)
- ▶ Process overview
 - Scorecard development and refinement
 - Selecting starting “candidate KBP model”
 - Comparing “candidate model” scores, “manual plan” scores, and “final model”
- ▶ Automation tools
 - PlanScoreCard batch scoring
 - Structure generation
- ▶ Future Developments

3

3

Knowledge Based Planning Concepts

4

4

What is Knowledge Based Planning?


RapidPlan is a KBP product

A method for predicting DVH's which utilizes Machine Learning

- Based on patient geometry
- Previous knowledge contained in a set of existing patient "training set" cases

Optimization objectives are automatically populated at the lower level of estimated DVH bands

What relative weight/priority should be used?



5

Objective Priorities


Setting the priorities for auto-created optimization objectives is not a straightforward process

Traditionally start with a "good guess" for auto-created optimization objectives for the model to act as "a starting point"

- With careful tuning of the optimization objectives, RapidPlan can create high quality plans -- with a single button press
- **Our method was to tune the auto-created optimization objectives via a ScoreCard**

6

DVH estimation is straightforward
 Consistently planned, high quality cases will help eliminate outliers




7

Head and Neck Cases are Heterogenous

Large variation in:


- The number of targets and prescription levels
 - Target geometry (size and shape)
 - Target location (proximity and overlap with OAR's)
- ▶ Variability requires different planning techniques and specific guidance to achieve desired plan quality
- ▶ Manual planning is time consuming and has a wide range of results based on the skill of the treatment planner




8

Head and Neck Cases are Heterogenous

Variation in high dose PTV size and location



Larynx




Oropharynx


9

Head and Neck Cases are Heterogenous

Variation in laterality, number of PTV's, and prescription dose



Unilateral: SIB 66, 60, & 56Gy



Bilateral: SIB 63, 60, 57 & 54Gy

10

10

Process

- **Theory**
 - Knowledge based planning models created with heterogenous cases can lead to wide range of DVH estimations
 - These wider DVH estimations can reduce plan quality

- **Goal**
 - Create a specific model for SIB (70, 63, & 56Gy) bilateral H&N with **both parotids** partially involved that requires minimal (if any) interaction from a dosimetrist

- **Method**
 - Test several existing models to select the best candidate for initial training set case generation

11


11

Head and Neck Model Background

12


12

AI Clustered Plans by Parotid Dose



13

Plan Clustering






Parotid Gland DVH

14

14

Data Collection Process

1. Import and evaluate the 118 cases that fall in C0 (dose to both parotids)
2. Keep cases that are bilateral with three targets (27 total)
3. Add any missing OAR contours and delineate the ipsilateral and contralateral parotids
4. Use an existing KBP model as a starting point for optimization
(clinical doses are not used)
5. Design scorecard for evaluation and C0 model tuning



15

15

Scorecard Development

16

16

Scorecards: Basic Principles

Between two different patients, scores should not be directly compared

Designed to analyze multiple plan options for the same patient/case

Scorecards can be created for a single patient, but better to be created per protocol or class of patients

Max score should not be achievable

- Scorecard total
- Per metric


Powerful tool for retrospective plan quality analysis

- Create dosimetric wish list to document best achievable today



H&N Scorecard Development

Protocols lack specificity



Toward Systematic Assessment and Improvement of Radiation Therapy Plan Quality of Cooperative Group Trial Submissions: A Report From the Children's Oncology Group

Arthur J. Olch, PhD,^{a,*} Mahesh Gopalakrishnan, MS,^b Erin S. Murphy, MD,^c Shannon M. MacDonald, MD,^d and Chia-Ho Hua, PhD^e

^aDepartment of Radiation Oncology, University of Southern California and Children's Hospital Los Angeles, Los Angeles, California; ^bDepartment of Radiation Oncology, Northwestern Memorial Hospital, Chicago, Illinois; ^cDepartment of Radiation Oncology, Cleveland Clinic, Cleveland, Ohio; ^dDepartment of Radiation Oncology, Massachusetts General Hospital, Boston, Massachusetts; and ^eDepartment of Radiation Oncology, St. Jude Children's Research Hospital, Memphis, Tennessee

Received 6 December 2022; accepted 2 February 2023

Purpose: This study evaluates the quality of plans used for the treatment of patients in the Children's Oncology Group study ACNS1123. Plan quality is quantified based on a scoring system specific to the protocol. In this way, the distribution of plan quality scores is determined that can be used to identify plan quality issues for this study and for future plan quality improvement.

Methods and Materials: ACNS1123 stratum 1 patients (70) were evaluated. This included 50 photon and 20 proton plans. Digital Imaging and Communications in Medicine (DICOM) structure and dose data were obtained from the Children's Oncology Group. A commercially available plan quality scoring algorithm was used to create a scoring system we designed using the protocol dosimetric requirements. The whole ventricle and boost planning target volumes (PTVs) could earn a maximum of 70 points, whereas the organs at risk could earn 30 points (total maximum score of 100 points). The scoring algorithm adjusted scores based on the difficulty in achieving the structure dose requirements, which depended on the proximity of the PTVs and the dose gradients achieved relative to the organs at risk. The distribution of plan scores was used to determine the mean, median, and range of scores.

Results: The median adjusted plan quality scores for the 20 proton and 50 photon plans were 83.3 and 86.9, respectively. The range of adjusted scores (maximum to minimum) was 50 points. The average score adjustment was 7.4 points. Photon and proton plans performed almost equally. Average plan quality by individual structure revealed that the brain stem, PTV boost, and cochlea lost the most points.

Conclusions: This report is the first to systematically analyze overall radiation therapy plan quality scores for an entire cohort of patients treated in a cooperative group clinical trial. The methodology demonstrated a large variation in plan quality in this trial. Future clinical trials could potentially use this method to reduce plan quality variability, which may improve outcomes.

© 2023 American Society for Radiation Oncology. Published by Elsevier Inc. All rights reserved.

H&N Scorecard Development

Protocols lack specificity

LBA 02
Radiotherapy with Durvalumab vs. Cetuximab in Patients with Locoregionally Advanced Head and Neck Cancer and a Contraindication to Cisplatin: Phase II Results of NRG-HN004

L.K. Mell^{1,2}, B. Torres-Saavedra³, S. Wong⁴, S. Chang⁵, J.A. Kish⁶, A.J. Minn III⁷, R. Jordan⁸, T. Liu⁹, M.T. Truong¹⁰, E. Wilmquist¹¹, T. Wise-Draeger¹², C.P. Rodriguez¹³, A. Musaddiq¹⁴, B.M. Readle¹⁵, C. Henson¹⁶, S. Narayan¹⁷, S.A. Spencer¹⁸, J. Harris⁹, S.S. Yom¹⁹

Show more
 + Add to Mendeley Share Cite
<https://doi.org/10.1016/j.jrobp.2022.09.003> Get rights and content »

Purpose/Objective(s)

The optimal treatment for patients with locoregionally advanced head and neck squamous cell carcinoma (HNSCC) and contraindication to cisplatin is uncertain. This trial (NCT03258554) tested the primary hypothesis that radiation therapy (RT) with concurrent and adjuvant durvalumab, a PD-L1 inhibitor, improves progression-free survival (PFS) compared to standard RT with cetuximab.

Organ at Risk	Objective	
	D _{max} (cGy)	Dose-Volume Limit
SpinalCord_05	Myelopathy V(5000 cGy) < 0.03 cc	V(3000 cGy) < 45% V(4000 cGy) < 10% Lhermitte's < 2000
Brainstem_03	Neuropathy V(5400 cGy) < 0.03 cc	Or D(2.7 cc) < 5500 cGy D(0.9 cc) < 6000 cGy Nausea < 3000 cGy
PartialBrain Chiasm	V(5400 cGy) < 0.03 cc	V(5000 cGy) ALARA
Pituitary / Hypothalamus		Mean < 4500 cGy
Lacrimal Glands		Mean < 3500 cGy
OpticNerve_R	V(6000 cGy) < 0.03 cc	
OpticNerve_L	V(6000 cGy) < 0.03 cc	
Retina	V(5000 cGy) < 0.03cc	Mean < 4500 cGy
Cochlea		Nausea V(4000 cGy) < 80% Hearing loss < 4500 cGy
MiddleEar		Mean < 4600 cGy
Parotid_R		Mean < 2600 cGy
Parotid_L		Mean < 2600 cGy
Parotid_Stem_R		
Parotid_Stem_L		
Submandibula_R		Mean < 3900 cGy
Submandibula_L		Mean < 3900 cGy
OralCavity		Mean < 3200 cGy
Larynx		V(3500 cGy) < 79% V(4500 cGy) < 45% V(5500 cGy) < 32% V(6500 cGy) < 22% Aspiration < 4100 cGy PEG Dependence < 5100 cGy
Thyroid		10-15cc Gland - mean < 1000 cGy 20 cc Gland - mean < 2500 cGy 25 cc Gland - mean < 4000 cGy
SPC/MPC/IPC		Aspirabon < 5400 cGy Structure < 5400 cGy PEG Dependence < 5100 cGy
IPC		V(4000 cGy) < 65% V(5000 cGy) < 47%

H&N Scorecard Development

Metric Editor

The screenshot shows the Metric Editor interface with the following components:

- Table 1:** A table with columns RANK, VALUE, SCORE, VARIATION, and COLOR.

RANK	VALUE	SCORE	VARIATION	COLOR
0	10	12	<input type="checkbox"/>	Dark Green
1	18	11.5	<input type="checkbox"/>	Light Green
2	26	8	<input checked="" type="checkbox"/>	Yellow
3	35	0	<input type="checkbox"/>	Orange
- Table 2:** A table with columns Scorecard Structure Id, Plan Structure, and Metric Type.


Scorecard Structure Id	Plan Structure	Metric Type
ParotidIps-PTV	ParotidIps-PTV	MeanDose
- Table 3:** A table with columns Metric Comment and Value.

Metric Comment	Value
ALARA (10-18Gy) /NRG-HN004 (26Gy)	4
- Graph:** A line graph titled "Mean Dose Score" showing Score vs Dose [Gy]. The x-axis ranges from 10 to 35 Gy, and the y-axis ranges from 0 to 12. The graph shows a piecewise linear function connecting points (10, 12), (18, 11.5), (26, 8), and (35, 0). The segment between 10 and 18 Gy is dark green, 18-26 Gy is light green, and 26-35 Gy is yellow.

1. Points assigned for mean dose in Gy
2. Metric points scorecard expressed in qualitative colors
 - Orange = 0-8 points
 - Dark Green = Full points
3. Points plotted as a piecewise linear function
 - Yellow symbolizes variation range
4. Comment box used for referencing protocol

H&N Scorecard Development

Contra and Ipsilateral Parotid Metrics




1. Parotid glands were designated as either ipsilateral or contralateral
2. More total points are awarded for contralateral parotid
3. Contralateral parotid has a steeper variation point slope
4. Ipsilateral parotid has a higher failing dose (0 points)

21

H&N Scorecard Development


Point Slopes



22

H&N Scorecard Development

Point Slopes



23

23

H&N Scorecard Development

Point Distribution

Balance Target metrics (up to 50%)

- Max/min dose
- % coverage
- Homogeneity
- Conformality

OAR metrics

- Volume at Dose
- Mean dose
- Dose at Volume (0.03cc)

Highest points to PTV coverage

Add multiple dose level rings

(PTV56, PTV63, & PTV70)

Max score should not be achievable

- Rarely: single metric (out of field)
- Never: total score card

Target metrics (37.5%)

ID	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV70OPT	Volume at 70Gy [%]	New HN Model	97.35 %	19.12	20.00	
3	PTV70OPT	Dose at 0.03CC [Gy]	New HN Model	74.33 Gy	8.91	10.00	
47	RingPTV70	Dose at 0.03CC [Gy]	New HN Model	71.33 Gy	4.53	5.00	

OAR metrics (62.5%)

34	Mandible-PTV	Volume at 60Gy [%]	New HN Model	11.26 %	1.80	2.00	
29	Lips	MeanDose [Gy]	New HN Model	16.05 Gy	6.19	7.00	
37	Esophagus	Dose at 0.03CC [Gy]	New HN Model	52.07 Gy	3.60	3.00	

24

24

H&N Scorecard Development

53 Metrics with 260 total points

Metric Id	Structure Id	Metric	Max Score
0	PTV70OPT	Volume at 70Gy [%]	20
1	PTV70OPT	Dose at 99.5% [Gy]	1.5
2	PTV70OPT	Dose at 0.03CC [Gy]	10
3	PTV63	Volume at 63Gy [%]	17
4	PTV63	Dose at 99.5% [Gy]	1.5
5	PTV63-PTV70	Volume at 66.15Gy [%]	8
6	PTV56	Volume at 56Gy [%]	15
7	PTV56	Dose at 99.5% [Gy]	1.5
8	PTV56-PTV63	Volume at 58.8Gy [%]	8
9	SpinalCord_05	Dose at 0.03CC [Gy]	6.5
10	SpinalCord_05	Volume at 40Gy [%]	2
11	SpinalCord_05	Volume at 30Gy [%]	2
12	Brainstem_03	Dose at 0.03CC [Gy]	4
13	Brain	Dose at 0.03CC [Gy]	2
14	Brain	Volume at 50Gy [CC]	3
15	Pituitary	MeanDose [Gy]	1
16	Chiasm	Dose at 0.03CC [Gy]	3
17	OpticNerve_L	Dose at 0.03CC [Gy]	3
18	OpticNerve_R	Dose at 0.03CC [Gy]	3
19	LacrimalGlands	MeanDose [Gy]	3
20	Cochlea_R	Volume at 40Gy [%]	3
21	Cochlea_L	Volume at 40Gy [%]	3
22	Lens_R	Dose at 0.03CC [Gy]	2.5
23	Lens_L	Dose at 0.03CC [Gy]	2.5
24	Eye_R	Dose at 0.03CC [Gy]	2
25	Eye_R	MeanDose [Gy]	2
26	Eye_L	Dose at 0.03CC [Gy]	2
27	Eye_L	MeanDose [Gy]	2


Metric Id	Structure Id	Metric	Max Score
28	Lips	MeanDose [Gy]	7
29	ParotdIps-PTV	MeanDose [Gy]	12
30	ParotdCon-PTV	MeanDose [Gy]	15
31	PharConst-PTV	MeanDose [Gy]	5
32	Mandible-PTV	Volume at 70Gy [%]	5
33	Mandible-PTV	Volume at 60Gy [%]	2
34	Mandible-PTV	Volume at 50Gy [%]	2
35	Esophagus	MeanDose [Gy]	4
36	Esophagus	Dose at 0.03CC [Gy]	3
37	OCavity-PTV	MeanDose [Gy]	6
38	OCavity-PTV	Dose at 0.03CC [Gy]	2
39	Larynx-PTV	MeanDose [Gy]	7
40	Thyroid-PTV	MeanDose [Gy]	2
41	BrachialPlexus_L	Dose at 0.1CC [Gy]	4
42	BrachialPlexus_R	Dose at 0.1CC [Gy]	4
43	SubmandL-PTV	MeanDose [Gy]	9.25
44	SubmandR-PTV	MeanDose [Gy]	9.25
45	TMJoint	Dose at 0.03CC [Gy]	2
46	RingPTV70	Dose at 0.03CC [Gy]	5
47	RingPTV63	Dose at 0.03CC [Gy]	5
48	RingPTV56	Dose at 0.03CC [Gy]	5
49	Posterior_Neck	Volume at 35Gy [%]	5
50	Trachea	MeanDose [Gy]	2.5
51	Lungs	Volume at 20Gy [CC]	2
52	Shoulders	MeanDose [Gy]	1

25

H&N Scorecard Development

Highlights

- PTV70_OPT (Cropped from the Brachial Plexus)
- Lower PTV's cropped from the higher PTV's (i.e. PTV63 – PTV70)
- OARs were evaluated by using substructures not overlapping with PTVs
- Unique ring structures were added for each PTV
- Shoulder and posterior neck structures were added




26

Model Training Process

27


27

Traditional Model Training Process




28


28



29



30



31

Testing candidate KBP models

Test existing KBP models

Patient	Course	Plan	<input type="checkbox"/>	Score (X/250)
002445	RC	Candidate 1	<input checked="" type="checkbox"/>	181.82
002445	RC	Candidate 4	<input checked="" type="checkbox"/>	144.00
002445	RC	Candidate 3	<input checked="" type="checkbox"/>	161.09
002445	RC	Candidate 2	<input checked="" type="checkbox"/>	138.51

Special thanks to Vanessa Magliari for providing Candidate 1 (top candidate model) used throughout the rest of this work!


31	ParotdCon-PTV	MeanDose [Gy]	Candidate 1	18.78 Gy	10.83	15.00
			Candidate 4	22.67 Gy	5.00	
			Candidate 3	20.29 Gy	8.57	
			Candidate 2	33.28 Gy	0.00	

IDEAL[15] GOOD[12] VARIATION[0]

7	PTV56	Volume at 56Gy [%]	Candidate 1	98.08 %	14.68	15.00
			Candidate 4	89.43 %	0.00	
			Candidate 3	95.12 %	13.09	
			Candidate 2	59.09 %	0.00	

SUB-OPTIMAL[0] VARIATION[13] GOOD[14.5] IDEAL[15]

32



33


Addressing Target Conformity

Gather feedback from original plans


- **Goal:** Improve Target conformity while still sparing parotid glands.
- **Method:** Create individual ring structures for each PTV
- **Process:** Include the ring structures in the new model and scorecard

34


34




35




36




37




38




39




40




41



42



43



44

Training Set Cases (Edited Top Candidate Cases)

- Record the total score of all optimized plans
- Calculate the average score (used as a benchmark to track the overall performance of future versions of the model)
- Take note of any unique patient anatomy that could affect the overall score for future reference
- Utilize these cases to train the first iteration of the new model


Patient	Notes	Plan Score
RP_HN_C05_006	Thymus	230.66
RP_HN_C05_016	Thymus	208.04
RP_HN_C05_021	Thymus	233.95
RP_HN_C05_131	Oropharynx	215.58
RP_HN_C05_196	Posterior LNs	209.78
RP_HN_C05_271	Oropharynx	232.01
RP_HN_C05_293	Oropharynx	231.87
RP_HN_C05_314	Oropharynx	231.87
RP_HN_C05_318	Oropharynx	221.25
RP_HN_C05_325	Oropharynx	226.71
RP_HN_C05_342	Oropharynx	236.25
RP_HN_C05_345	Oropharynx	226.1
RP_HN_C05_347	Oropharynx	233.55
RP_HN_C05_356	Oropharynx	237.51
RP_HN_C05_402	Oropharynx	217.5
RP_HN_C05_413	Oropharynx	228.93
RP_HN_C05_434	Posterior LNs	232.1
RP_HN_C05_452	Inferior PTV70	231.82
RP_HN_C05_529	Anterior LNs	223.59
RP_HN_C05_604	OVL with mandible	215.48
RP_HN_C05_623	Oropharynx	220.29
RP_HN_C05_681	Oropharynx	229.21
RP_HN_C05_682	Oropharynx	230.86
RP_HN_C05_768	Small PTV63	234.14
RP_HN_C05_775	High OVL with BP	224
RP_HN_C05_785	Oropharynx	232.33
RP_HN_C05_771	Posterior LNs	232.59
Average Score		226.96

Score plans

45

45

Plan Score Comparison Top candidate model vs edited top candidate model



Score plans


Average Scores

Top Candidate Model – 210.15

Edited Top Candidate Model – 226.96

46

46




47

Create initial model from plans


Target ID	Vol (cc)	Dose	Priority	gt/UD/A
Line (preferring target)	Generated	Generated	75	X
Parallels-PTV (Parallels)				
Mean		Generated	65	X
Line (preferring target)	Generated	Generated	70	X
Parallels-PTV (Parallels)				
Mean		Generated	65	X
Line (preferring target)	Generated	Generated	60	X
Parallels (Parallels, 59798, 59797)				
Mean		Generated	100	X
Line (preferring target)	Generated	Generated	90	X
Parallels-PTV (PRV)				
Upper (fixed dose, generated vol.)	Generated	19,000 Gy	50	X

#	Patient ID/Case ID/Plan ID	Plan Description	Structure Matching	Included	Extracted	In Model
1	RP_HN_C06_006C1KBP006	70,000 Gy	Target: 33 Other: 4244	Yes	Yes	18.0.0 X
2	RP_HN_C06_016C1KBP016	70,000 Gy	Target: 33 Other: 4244	Yes	Yes	18.0.0 X
4	RP_HN_C06_131C1KBP131	70,000 Gy	Target: 33 Other: 4344	Yes	Yes	18.0.0 X
6	RP_HN_C06_196C1KBP196	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X
7	RP_HN_C06_271C1KBP271	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X
8	RP_HN_C06_293C1KBP293	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X
9	RP_HN_C06_314C1KBP314	70,000 Gy	Target: 33 Other: 4344	Yes	Yes	18.0.0 X
10	RP_HN_C06_318C1KBP318	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X
11	RP_HN_C06_325C1KBP325	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X
12	RP_HN_C06_342C1KBP342	70,000 Gy	Target: 33 Other: 4444	Yes	Yes	18.0.0 X


48




49




50




51



52



53




54

Initial Results

AAMD 2023 Plan Study- Isodose


- “Oral avoid” structure was used for plan challenge case (shown in pink)
- Added “oral cavity” structure (purple) that matched the model training set cases.
- Unwanted dose bridging still occurred within the “oral avoid”
- Manual adjustments were needed during optimization to pass



55

Initial Results


AAMD 2023 Plan Study- DVH



56

Initial Results

AAMD 2023 Plan Study- Scorecard (133.03 / 150 points)



57

Plan Score Comparison

Total Scores: Top candidate vs current model (Four external validation cases)




Average Scores
 Top Candidate model – 205.4
 Current Model – 222.06

58

Initial Results


Improved Conformality

Ring PTV56




Cropped 0.2cm from PTV56, 0.6cm from PTV63, and 0.9cm from PTV70

Ring PTV63



Cropped 0.2cm from PTV63 and 0.6cm from PTV70

Ring PTV70




Cropped 0.2cm from PTV70

59

Initial Results

Improved parotid sparing and DVHe



60


Automation Tools

61

61

Structure Generation

Automatically create rings



BUILD A STRUCTURE

AND OR SUB Edit Groupings

Operation	Structure	Margin [mm]	
	PTV56	- 30	+
SUB	PTV56	- 2	+
SUB	PTV63	- 6	+
SUB	PTV70	- 9	+
AND	BODY	- 0	+

Group1
 <PTV56>|30 SUB <PTV56>|2 SUB <PTV63>|6 SUB <PTV70>|9 AND <BO

Group1
 - 0 +

Structure Builder Comment
 [<PTV56>|30 SUB <PTV56>|2 SUB <PTV63>|6 SUB <PTV70>|9 AND <BODY>]

Structure Id: RingPTV56 FINALIZE STRUCTURE

62

62

Structure Generation

Boolean operations to generate rings are nested in the scorecard metrics

```

"Structure": {
  "StructureId": "RingPTV56",
  "StructureComment": "({<PTV56>|30| SUB <PTV56>|2 SUB <PTV63>|6 SUB <PTV70>|9 AND <BODY>})",
  "TemplateStructureId": "RingPTV56",
  "StructureCode": "9999",
  "AutoGenerated": true,
  "bMakeVisibleInPatientSearch": false
},
"Rank": 0,
"MetricType": "DoseAtVolume",
"InputValue": 0.03,
"InputUnit": "CC",
"OutputUnit": "Gy",
"AcceptableLevel": null,
"VariationLevel": null,
"FailLevel": null,
"HI_HiValue": 0.0,
"HI_LowValue": 0.0,
"HI_Target": 0.0,
"HI_TargetUnit": null,
"MetricComment": null,
"TemplateNumber": 48,
"ScorePoints": [
  {
    "PointX": 56.0,
    "Score": 5.0,
    "Variation": false,
    "Colors": [
      0.0,
      153.0,
      0.0
    ],
    "Label": "IDEAL[5]",
    "ColorValue": 5.0
  }
]

```

63

63

Batch Scoring

Score multiple plans at once

PLAN SCORECARD (NOT VALIDATED FOR CLINICAL USE)

SCORE CARD SELECTION

Score Card ID: Plan R matches scorecard.
H&N-KBP-Draft0.json

PLAY SELECTION

Patient ID: RP_HN_C05_131

Patient	Course	Plan	Score (X/25C)
RP_HN_C05_131	C1	KBP131D2	221.12

Plan Scores: RP_HN_C05_131-[C1] KBP131D2: 221.12/250.75 (88.19%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV700PT	Volume at 700g [%]	KBP131D2	99.73 %	99.01	100.00	
2	PTV700PT	Dose at 99.5% [Gy]	KBP131D2	70.27 Gy	1.50	1.50	
3	PTV700PT	Dose at 0.03CC [Gy]	KBP131D2	74.09 Gy	9.08	95.00	
4	PTV63	Volume at 930g [%]	KBP131D2	96.61 %	16.65	17.00	
5	PTV63	Dose at 99.5% [Gy]	KBP131D2	62.13 Gy	1.33	1.50	
6	PTV63-PTV70	Volume at 66.150g [%]	KBP131D2	19.07 %	7.55	8.00	
7	PTV56	Volume at 560g [%]	KBP131D2	98.64 %	14.77	15.00	

Norm Export Print

Bounded by the terms of the [Vendor EULA](#)

64

Batch Scoring

Query TPS, select patients + course/plan, and match structures to scorecard

PATIENTSELECTIONVIEW

Patients

RP_HN_C05_785

RP_HN_C05_785, Cluster0 (3 Targets 54.4/64.6/68Gy) (RP)

Selected Patients

- RP_HN_C05_016
- RP_HN_C05_016, Cluster0 (3 Targets 54/60/66Gy)
- RP_HN_C05_785
- RP_HN_C05_785, Cluster0 (3 Targets 54.4/64.6/68Gy)
- RP_HN_C05_413
- RP_HN_C05_413, Cluster0 (3 Targets 48/54/60Gy)
- RP_HN_C05_342
- RP_HN_C05_342, Cluster0 (3 Targets 48/54/60Gy)
- RP_HN_C05_006
- RP_HN_C05_006, Cluster0 (3 Targets 54/60/66Gy)
- RP_HN_C05_021
- RP_HN_C05_021, Cluster0 (3 Targets 54/60/66Gy)
- RP_HN_C05_131
- RP_HN_C05_131, Cluster0 (3 Targets 54/60/66Gy)
- RP_HN_C05_196
- RP_HN_C05_196, Cluster0 (3 Targets 54.4/61.2/68Gy)
- RP_HN_C05_271
- RP_HN_C05_271, Cluster0 (3 Targets 48/54/60Gy)

Plan selection for RP_HN_C05_785

Scorecard imported with 59 metrics

- C1: CopyJIC.0
- C1: Hal4Arc
- C1: KBP785
- C1: KBP785D2

Validation for KBP016

Scorecard Structure	Matched Structure	OK	Dict Match	Local Match	Fix
	eye_n				Fix
Eye_L	Eye_L	✓			Fix
Lips	Lips	✓			Fix
ParotdIps-PTV	ParotdIps-PTV	✓			Fix
ParotdCon-PTV	ParotdCon-PTV	✓			Fix
PharConst-PTV	PharConst-PTV	✓			Fix
Mandible-PTV	Mandible-PTV	✓			Fix
Esophagus	Esophagus	✓			Fix
OCavity-PTV	OCavity-PTV	✓			Fix

SCORE SELECTIONS
SAVE PATIENT SELECTIONS
REMOVE SELECTED PATIENT
CANCEL SELECTIONS

65

Batch Scoring

Contralateral parotid scores for 27 plans

Contralateral Parotid (Mean Dose)

- 27 plans scored
- Value = Mean dose
- Score out of 15 possible points

Plan Id	Value	Score	Max	Metric Plot
KBP016D2	10.90 Gy	13.64	15.00	
KBP785D2	11.63 Gy	13.47	15.00	
KBP413D2	12.89 Gy	13.18	15.00	
KBP342D2	7.95 Gy	14.32	15.00	
KBP006D2	5.13 Gy	14.97	15.00	
KBP021D2	8.37 Gy	14.22	15.00	
KBP131D2	19.15 Gy	10.27	15.00	
KBP196D2	15.45 Gy	12.59	15.00	
KBP271D2	8.18 Gy	14.27	15.00	
KBP293D2	14.21 Gy	12.87	15.00	
KBP314D2	10.32 Gy	13.77	15.00	
KBP318D2	12.52 Gy	13.27	15.00	
KBP325D2	13.28 Gy	13.09	15.00	
KBP345D2	17.63 Gy	12.09	15.00	
KBP347D2	8.69 Gy	14.15	15.00	
KBP356D2	9.37 Gy	13.99	15.00	
KBP402D2	10.62 Gy	13.70	15.00	
KBP434D2	9.53 Gy	13.96	15.00	
KBP452D2	11.28 Gy	13.55	15.00	
KBP529D2	16.09 Gy	12.44	15.00	
KBP604D2	19.50 Gy	9.75	15.00	
KBP623D2	12.52 Gy	13.26	15.00	
KBP681D2	14.28 Gy	12.86	15.00	
KBP682D2	12.82 Gy	13.20	15.00	
KBP768D2	15.92 Gy	12.48	15.00	

Metric Plot

- Max points: 5Gy
- Variation: 18Gy
- Zero points 26Gy

66

Batch Scoring

Score results can be exported as a CSV for data analysis

PLAN SCORECARD ***NOT VALIDATED FOR CLINICAL USE***

SCORE CARD SELECTION
 Score Card ID: Plan Ra matches scorecard.
 HSN-KBP-Draft02.json

PLAN SELECTION
 Patient ID: EP_HN_C05_001

Patient	Course	Plan	Score (X/256)
RP_HN_C05_293	C1	KBP293D2	234.63
RP_HN_C05_314	C1	KBP314D2	231.05
RP_HN_C05_318	C1	KBP318D2	224.58
RP_HN_C05_325	C1	KBP325D2	220.57
RP_HN_C05_345	C1	KBP345D2	223.49
RP_HN_C05_347	C1	KBP347D2	236.09
RP_HN_C05_356	C1	KBP356D2	237.87
RP_HN_C05_402	C1	KBP402D2	222.80
RP_HN_C05_434	C1	KBP434D2	234.95
RP_HN_C05_452	C1	KBP452D2	232.84
RP_HN_C05_529	C1	KBP529D2	225.67
RP_HN_C05_604	C1	KBP604D2	225.69
RP_HN_C05_623	C1	KBP623D2	222.99
RP_HN_C05_681	C1	KBP681D2	232.68
RP_HN_C05_682	C1	KBP682D2	232.58
RP_HN_C05_768	C1	KBP768D2	236.20
RP_HN_C05_771	C1	KBP771D2	233.73
RP_HN_C05_775	C1	KBP775D2	229.40

Plan Scores: RP_HN_C05_016-[C1] KBP016D2: 212.66/258.00 (82.43%)
 RP_HN_C05_785-[C1] KBP785D2: 235.42/258.00 (91.25%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV700PT	Volume at 70Gy [L]	KBP016D2	97.69 %	19.23	20.00	
			KBP785D2	88.52 %	19.51	Score Stats: Mean=19.84, Minus=17.92	
			KBP413D2	98.71 %	19.57		
			KBP342D2	98.52 %	19.51		
			KBP006D2	99.56 %	19.85		
			KBP021D2	99.37 %	19.79		
			KBP131D2	99.73 %	19.91		
			KBP196D2	98.33 %	19.44		
			KBP271D2	98.92 %	19.64		
			KBP293D2	98.36 %	19.79		
			KBP314D2	98.55 %	19.52		
			KBP318D2	97.22 %	19.07		
			KBP325D2	99.42 %	19.81		
			KBP345D2	99.05 %	19.68		
			KBP347D2	95.92 %	17.92		
			KBP356D2	98.96 %	19.65		
			KBP402D2	97.91 %	19.30		
			KBP434D2	99.42 %	19.81		
			KBP452D2	98.68 %	19.56		
			KBP529D2	99.09 %	19.70		
			KBP604D2	97.70 %	19.23		
			KBP623D2	99.82 %	19.94		
			KBP681D2	97.81 %	19.27		
			KBP682D2	97.82 %	19.27		
			KBP768D2	88.73 %	19.58		

Buttons: Norms, Export, Print

67

Batch Scoring

Contralateral parotid scores for 27 plans in CSV

Metric Id	Patient Id	Course Id	Plan Id	StructureId	Metric Text	Metric Value	Score	Max Score
30	RP_HN_C05_016	C1	KBP016D2	ParotdCon-PTV	MeanDose [Gy]	10.89995406	13.63847214	15
30	RP_HN_C05_785	C1	KBP785D2	ParotdCon-PTV	MeanDose [Gy]	11.63374166	13.46913654	15
30	RP_HN_C05_413	C1	KBP413D2	ParotdCon-PTV	MeanDose [Gy]	12.88846921	13.17958403	15
30	RP_HN_C05_342	C1	KBP342D2	ParotdCon-PTV	MeanDose [Gy]	7.948892451	14.31948636	15
30	RP_HN_C05_006	C1	KBP006D2	ParotdCon-PTV	MeanDose [Gy]	5.12981632	14.97004239	15
30	RP_HN_C05_021	C1	KBP021D2	ParotdCon-PTV	MeanDose [Gy]	8.373056649	14.22160231	15
30	RP_HN_C05_131	C1	KBP131D2	ParotdCon-PTV	MeanDose [Gy]	19.15231413	10.27152888	15
30	RP_HN_C05_196	C1	KBP196D2	ParotdCon-PTV	MeanDose [Gy]	15.44828304	12.58885776	15
30	RP_HN_C05_271	C1	KBP271D2	ParotdCon-PTV	MeanDose [Gy]	8.182532946	14.26556932	15
30	RP_HN_C05_293	C1	KBP293D2	ParotdCon-PTV	MeanDose [Gy]	14.21413461	12.87366124	15
30	RP_HN_C05_314	C1	KBP314D2	ParotdCon-PTV	MeanDose [Gy]	10.31891416	13.77255827	15
30	RP_HN_C05_318	C1	KBP318D2	ParotdCon-PTV	MeanDose [Gy]	12.51621788	13.26548818	15
30	RP_HN_C05_325	C1	KBP325D2	ParotdCon-PTV	MeanDose [Gy]	13.27648483	13.09004196	15
30	RP_HN_C05_345	C1	KBP345D2	ParotdCon-PTV	MeanDose [Gy]	17.6276355	12.08593027	15
30	RP_HN_C05_347	C1	KBP347D2	ParotdCon-PTV	MeanDose [Gy]	8.690954828	14.14824119	15
30	RP_HN_C05_356	C1	KBP356D2	ParotdCon-PTV	MeanDose [Gy]	9.367021136	13.99222589	15
30	RP_HN_C05_402	C1	KBP402D2	ParotdCon-PTV	MeanDose [Gy]	10.62072982	13.7029085	15
30	RP_HN_C05_434	C1	KBP434D2	ParotdCon-PTV	MeanDose [Gy]	9.527350758	13.95522675	15
30	RP_HN_C05_452	C1	KBP452D2	ParotdCon-PTV	MeanDose [Gy]	11.2752259	13.55187095	15
30	RP_HN_C05_529	C1	KBP529D2	ParotdCon-PTV	MeanDose [Gy]	16.08590426	12.4417144	15
30	RP_HN_C05_604	C1	KBP604D2	ParotdCon-PTV	MeanDose [Gy]	19.49747593	9.753786099	15
30	RP_HN_C05_623	C1	KBP623D2	ParotdCon-PTV	MeanDose [Gy]	12.52315812	13.26388659	15
30	RP_HN_C05_681	C1	KBP681D2	ParotdCon-PTV	MeanDose [Gy]	14.28048378	12.8583499	15
30	RP_HN_C05_682	C1	KBP682D2	ParotdCon-PTV	MeanDose [Gy]	12.81607352	13.19629073	15
30	RP_HN_C05_768	C1	KBP768D2	ParotdCon-PTV	MeanDose [Gy]	15.91807267	12.48044477	15
30	RP_HN_C05_771	C1	KBP771D2	ParotdCon-PTV	MeanDose [Gy]	10.3761298	13.75935466	15
30	RP_HN_C05_775	C1	KBP775D2	ParotdCon-PTV	MeanDose [Gy]	19.93124133	9.103138009	15
Average						12.53779	13.0452	15

Head and neck batch scoring:

- 53 metrics per patient
- 27 patients
- 3 iterations of plan scoring

4293 rows of data

68

Normalize Plans to Max Score

PLAN SCORECARD (NOT VALIDATED FOR CLINICAL USE)

SCORE CARD SELECTION
 Score Card ID: Plan Rx matches scorecard.
 h5n-KBP-Draft01.json

PATIENT SELECTION
 Patient ID: RP_HN_CS_131
 Patient: RP_HN_CS_131 Course: C1 Plan: KBP131D2 Score D25C

Plan Scores: RP_HN_CS_131: [C1] KBP131D2: 221.12/250.75 (88.19%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV700PT	Volume at 70Gy [%]	KBP131D2	99.73 %	19.91	20.00	
2	PTV700PT	Dose at 99.5% [Gy]	KBP131D2	70.27 Gy	1.50	1.50	
3	PTV700PT	Dose at 0.03CC [Gy]	KBP131D2	74.09 Gy	9.08	10.00	
4	PTV63	Volume at 63Gy [%]	KBP131D2	98.61 %	16.65	17.00	
5	PTV63	Dose at 99.5% [Gy]	KBP131D2	62.13 Gy	1.33	1.50	
6	PTV63-PTV70	Volume at 66.15Gy [%]	KBP131D2	19.07 %	7.55	8.00	
7	PTV56	Volume at 56Gy [%]	KBP131D2	95.64 %	14.77	15.00	

Norm Export Print

69

Normalize Plans to Max Score

OPTIMIZING ON PLAN NORMALIZATION

Score at 116 = -5819.86
 Score at 118 = -5819.17
 Score at 98 = 202.11
 Score at 98.2 = 204.8
 Score at 98.4 = 207.18
 Score at 98.6 = 209.46
 Score at 98.8 = 211.71
 Score at 99 = 213.92
 Score at 99.2 = 215.69
 Score at 99.4 = 217.4
 Score at 99.6 = 218.95
 Score at 99.8 = 220.09
 Score at 100 = 221.11
 Score at 100.2 = 222.02
 Score at 100.4 = 222.77
 Score at 100.6 = 223.21
 Score at 100.8 = 222.58
 Score at 101 = 221.32
 Score at 101.2 = 209.51
 Score at 101.4 = -1816.67
 Score at 101.6 = -2825.88
 Score at 101.8 = -2825.63
 Score at 102 = -2825.38
 Score at 100.4 = 222.77
 Score at 100.42 = 222.83
 Score at 100.44 = 222.92
 Score at 100.46 = 222.98
 Score at 100.48 = 223.02
 Score at 100.5 = 223.09
 Score at 100.52 = 223.08
 Score at 100.54 = 223.12
 Score at 100.56 = 223.15
 Score at 100.58 = 223.18
 Score at 100.6 = 223.21
 Score at 100.62 = 223.23
 Score at 100.64 = 223.24
 Score at 100.66 = 223.25
 Score at 100.68 = 223.25
 Score at 100.7 = 223.25
 Score at 100.72 = 223.26
 Score at 100.74 = 223.22
 Score at 100.76 = 223.16
 Score at 100.78 = 223.08
 Score at 100.8 = 222.98

Max Score 223.257

ScoreCard Normalization: 100.72

* Activate Plan *
 CourseID: N-Opt
 PlanID: KBP131D2

Thousands of normalization iterations are used to find the max score.

Optimized Normalization: 100.72%

New Score: 223.26 points

70

Normalize Plans to Max Score

Plan copied under new course with Max

PLAN SCORECARD ***NOT VALORIED FOR CLINICAL USE***

SCORE CARD SELECTION

Score Card ID: Plan Rx matches scorecard.
H&N-KBP-Dual9 (jion)

Patient ID: RP_HN_C05_131

Patient	Course	Plan	Score DV25C
RP_HN_C05_131	CT	KBP131D2	221.12
RP_HN_C05_131	N-Gal	KBP131D2	223.6

Plan Scores: RP_HN_C05_131: [N-Gal] KBP131D2: 223.26/250.75 (89.84%)
RP_HN_C05_131: [CT] KBP131D2: 221.12/250.75 (88.19%)

Id	Structure	Score Metric	Plan Id	Value	Score	Max	Metric Plot
1	PTV700PT	Volume at 70Gy [%]	KBP131D2	98.92 %	19.64	20.00	
			KBP131D2	99.73 %	19.91	20.00	
SUB-OPTIMAL[0] VARIATION[14.57] GOOD[16.25] IDEAL[17]							
2	PTV700PT	Dose at 99.5% [Gy]	KBP131D2	66.77 Gy	1.46	1.50	
			KBP131D2	70.27 Gy	1.50	1.50	
SUB-OPTIMAL[0] VARIATION[0] GOOD[1.46] IDEAL[1.5]							
3	PTV700PT	Dose at 0.03CC [Gy]	KBP131D2	73.56 Gy	9.46	10.00	
			KBP131D2	74.09 Gy	9.08	10.00	
SUB-OPTIMAL[0] VARIATION[0] GOOD[9.08] IDEAL[10]							
4	PTV63	Volume at 63Gy [%]	KBP131D2	97.02 %	16.26	17.00	
			KBP131D2	98.61 %	16.65	17.00	
SUB-OPTIMAL[0] VARIATION[14.57] GOOD[16.25] IDEAL[17]							
5	PTV63	Dose at 99.5% [Gy]	KBP131D2	61.69 Gy	1.24	1.50	
			KBP131D2	62.13 Gy	1.33	1.50	
SUB-OPTIMAL[0] VARIATION[0] GOOD[1.24] IDEAL[1.5]							

Norm Report Print

Bound by the terms of the [VariantUSA](#)

71

Normalize Plans to Max Score

Over 2 points gained (221.12 vs 223.6)

Bold Crosshair (Normalized Plan)

- Lower PTV63 Coverage
- Improved PTV 56 heterogeneity
- Improved PTV 56 conformity

4	PTV63	Volume at 63Gy [%]	KBP131D2	97.02 %	16.26	17.00	
			KBP131D2	98.61 %	16.65	17.00	
SUB-OPTIMAL[0] VARIATION[14.57] GOOD[16.25] IDEAL[17]							
9	PTV56-PTV63	Volume at 58.8Gy [%]	KBP131D2	30.25 %	7.48	8.00	
			KBP131D2	42.08 %	6.59	8.00	
IDEAL[8] GOOD[7.5] GOOD[6] VARIATION[0]							
48	RingPTV56	Dose at 0.03CC [Gy]	KBP131D2	56.98 Gy	4.56	5.00	
			KBP131D2	57.39 Gy	3.78	5.00	
IDEAL[5] GOOD[4.3] VARIATION[0]							

Score Stats
Max=16.65
Mean=16.45
Min=16.26

Score Stats
Max=7.48
Mean=7.04
Min=6.59

Score Stats
Max=4.56
Mean=4.17
Min=3.78

72

Future Developments: Works in Progress


73

73

Additional Models

Unilateral only model with two PTV targets

Works in Progress




74

Additional Models

One parotid completely involved

Works in Progress



75

3D Dose Prediction Research

Works in Progress

- ▶ Predict 3D dose maps from individual patient data
- ▶ Use a CNN model that combines the 3D patient images with the used field geometry
- ▶ Input data: CT, PTVs, OARs, field geometry, and prescription dose levels






Figure 2. The schematic architecture of the DCNN, which consists of an encoder block, a dense feature aggregation block, and a decoder block. The numbers in the middle block represent the dilation rate of each convolution operation. The input is a four-channel feature image, and the output is the predicted dose image.

"Predicting Voxel-level Dose Distributions for Esophageal Radiotherapy Using Densely Connected Network with Dilated Convolutions", Zhang et al., Phys. Med. Biol. July 2020 (Anhui Univ.)

76



77



78

3D Dose Prediction Comparator for scorecard evaluation

Works in Progress

Scoring 3D generated doses quickly helps compare the performance of KBP models



79

Score based plan reoptimization Assign RapidPlan model to multiple cases

Works in Progress

1 Preview Plans

- 2019RP-EsoPh-Prelim5 - C1 - Plan2
- Test Case 1 - C1 - E6-3Haze4
- Test Case 2 - C1 - EsoPhoVnlan
- Test Case 3 - C1 - Plan17

2 Normal Tissue Objective Settings

- Auto NTO: Priority: 100
- Manual NTO: Priority: 100
- Distance from Surface [mm]: 0
- Start Percentage: 0
- End Percentage: 0
- Fall-off: 0

3 RapidPlan Model Settings

Target	RP Structure Id	Matched Structures	RP Target	Target Dose
<input checked="" type="checkbox"/>	PTV	PTV	PTV	0
<input checked="" type="checkbox"/>	PTVnoHeart	PTVnoHeart	PTVnoHeart	0
<input checked="" type="checkbox"/>	Heart	HeartHEARTs		
<input type="checkbox"/>	Kidneys	Kidney_L;Kidney_R;KIDNEY_L;KIDNEY_R		
<input type="checkbox"/>	Liver	Liver		
<input type="checkbox"/>	Lungs-PTV	Lungs-PTV;LUNG TOTAL-PTV		
<input type="checkbox"/>	PTV03Ring	PTV03Ring		
<input type="checkbox"/>	SpinalCanal	SPINAL CORD		

1. List of plans automatically selected from batch scoring template.
2. Specific NTO settings can be set. Additional optimization settings such as convergence mode, aperture shape control, and intermediate dose calculations will also be added.
3. Select the RapidPlan model and match the structures to the model.

80

Score based plan reoptimization

Run RapidPlan optimization on multiple cases

Works in Progress

RAPID PLANIMATOR

Progress...

Calculating DVH Estimates based on Esophagus2019PlanChallengeFinal

Score for 2019RP-Esoph-Prelim5 Course = C1 Plan = Plan2				Score for 2019RP-Esoph-Prelim5 Course = AutoCourse Plan = Plan3			
Score	Score Details	Variations	Flagged	Score	Score Details	Variations	Flagged
97.00 (64.67%)		0	6	93.06 (62.04%)		0	6

Score for Test Case 1 Course = C1 Plan = E6-3Haces4				Score for Test Case 1 Course = AutoCourse Plan = E6-2Haces4			
Score	Score Details	Variations	Flagged	Score	Score Details	Variations	Flagged
134.48 (89.65%)		1	0	134.48 (89.65%)		1	0

Score for Test Case 2 Course = C1 Plan = EsophageoVarian			
Score	Score Details	Variations	Flagged
71.04 (47.36%)		1	7

Score for Course = C1 Plan = Plan17			
Score	Score Details	Variations	Flagged

1. Calculates the DVHe from the selected model
2. Automatically proceeds with optimization, calculation, and plan scoring before starting the next case in the list
3. The newly optimized plans are then created under a new course
4. Plans can then be scored to evaluate any new changes that were made to a RapidPlan model
5. This process saves time when iteratively tuning models

81

81

Resources

82

82

Resources

RapidPlan Model

List of publicly shared RapidPlan Models, eventually including the final head and neck model discussed in this presentation

<https://medicalaffairs.varian.com/rapidplan-dir>

Publicly shared RapidPlan models

Please feel free to share your model via PeerAnswers and we will add it the next time this list is updated

Tx Site	Institution	PDF Publication	Rx (max/fom)/Fraction	Protocol Link	Share Link	Modality	Version	
Brain-HSWBv2	Varian MedAffairs	yes	30Gy	NRG-C001	VarianMedAffairs	Photon	15.6 + 17.0	
Brain-HSWBv1	Varian MedAffairs	yes	30Gy	NRG-C001	VarianMedAffairs	Photon	15.5	
Brain-MishHSWB200y	Varian MedAffairs	yes	20Gy	40y	N/A	TBD	15.6 + 17.0	
Brain-MoshHSWB20-800y	Varian MedAffairs	yes	20Gy	40y	CCTO CE.7	TBD	15.6 + 18.0	
Breast-LN	RHS-NorthernSydneyCC	yes	57Gy (5B)	2.28Gy	NCT02284733	VarianMedAffairs	Photon	16.1
Prostate	University of Miami	no	80Gy (5B)	20y	similar Rx released	Peer Answer	Photon	13.2,23
R Lung SBRT	UC San Diego	yes	56Gy	8Gy	similar Rx released	Peer Answer	Photon	13.2,23
Head-and-Neck	UC San Diego	yes	70Gy	20y	RT00819	Peer Answer	Photon	13.2,23
Prostate Fossa	UC San Diego	yes	70.20y (5-B)	1.80y	RT05P012	Peer Answer	Photon	13.2,23
L Lung SBRT	UC San Diego	yes	56Gy	8Gy	similar Rx released	Peer Answer	Photon	13.2,23
Prostate SV	UC San Diego	yes	81Gy	1.80y	similar Rx released	VarianMedAffairs	Photon	13.2,23
Prostate SV LN	University of Miami	no	80Gy (5B)	20y	similar Rx released	VarianMedAffairs	Photon	13.2,23
Prostate SV LN	University of Miami	no	80Gy (5B)	20y	similar Rx released	VarianMedAffairs	Photon	13.2,23
Esophagus	Varian MedAffairs	yes	50.40y	1.80y	RT00819	Peer Answer	Photon	16.1
Esophagus	Varian MedAffairs	yes	50.40y	1.80y	see publication	VarianMedAffairs	Photon	16.1
Liver	Varian MedAffairs	yes	50.40y	1.80y	NRG-0102	VarianMedAffairs	Photon	16.1
Head-and-Neck	Varian MedAffairs	no	70Gy (5B)	2.10y	RT00819	VarianMedAffairs	Photon	16.1
Head-and-Neck	Northampton General	no	65Gy (5B)	2.20y	RT00824	Peer Answer	Photon	15.0,24
Cervical cancer	Royal Surrey County	yes	50.40y	1.80y	INTERTEC_INTERFACE	Peer Answer	Photon	13.1,15
Cervix node inv	Tata Memorial Centre	no	55Gy (5B)	2.20y	EMBRACE II	VarianMedAffairs	Photon	13.1,15
Cervix node inv	Tata Memorial Centre	no	55Gy (5B)	2.20y	EMBRACE II	VarianMedAffairs	Photon	13.1,15
Head-and-Neck	Univ College London	no	55Gy (5B)	20y	RT00824	Peer Answer	Photon	15.0,24
Anorectal	UC San Diego/OribRT	yes	62Gy (5B)	2.20y	APROOTE	OribRT	Photon	15.6
Breast APBI	UC San Diego/OribRT	no	30Gy	30y	MSRP-BB RT00041	OribRT	Photon	16.1
Oligodendroma	UC San Diego/OribRT	no	60Gy (5B)	20y	RT00826	OribRT	Photon	15.6
RTN Multi-target	UC San Diego/OribRT	no	50.40y (5B)	1.80y	RT00818	OribRT	Photon	15.6
OVN single target	UC San Diego/OribRT	no	50.40y	1.80y	RT00818	OribRT	Photon	15.6
Head-and-Neck	UC San Diego/OribRT	no	70Gy	20y	RT00815	OribRT	Photon	16.1
Liver SBRT	UC San Diego/OribRT	no	60Gy	20Gy	similar Rx released	OribRT	Photon	15.6
Lung/Mediastinum	UC San Diego/OribRT	no	60Gy	20y	RT00817	OribRT	Photon	15.6
Lung SBRT (left)	UC San Diego/OribRT	yes	56Gy	8Gy	similar Rx released	OribRT	Photon	15.6
Lung SBRT (right)	UC San Diego/OribRT	yes	56Gy	8Gy	similar Rx released	OribRT	Photon	15.6
Lymphoma	UC San Diego/OribRT	no	30.60y	1.80y	similar Rx released	OribRT	Photon	15.6
Oligomet SBRT	UC San Diego/OribRT	no	45Gy	10Gy	NRG-0101	OribRT	Photon	15.6
Pancreas SBRT	UC San Diego/OribRT	no	50Gy	10Gy	NCT021044	OribRT	Photon	15.6
Prostate Fossa	UC San Diego/OribRT	yes	70.20y (5-B)	1.80y	similar Rx released	OribRT	Photon	15.6
Prostate Fossa	UC San Diego/OribRT	yes	70.20y	1.80y	RT00P12	OribRT	Photon	15.6
Prostate + Nodes	UC San Diego/OribRT	no	60Gy (5B)	30y	CHM	OribRT	Photon	15.6
Sarcoma	UC San Diego/OribRT	no	50Gy	20y	RT00030	OribRT	Photon	15.6
Spine SBRT	UC San Diego/OribRT	no	30Gy	6Gy	similar Rx released	OribRT	Photon	15.6

83

Resources

Scorecards

List of publicly shared scorecards, eventually including the final head and neck scorecard shown in this presentation

<https://medicalaffairs.varian.com/dose-scorecards>

Example Dosimetric ScoreCards

These example dosimetric scorecards are compatible with the **MAAS-PlanScoreCard** tool and are to be used for illustrative purposes only. Reference Protocols often provide only the clinical goal metrics, these example ScoreCards award points for exceeding base clinical goal thresholds. Varian does not provide medical advice, example ScoreCards are not medical advice.

Tx Site	Rx (Max/Non)	Fraction	Reference Protocol(s) / Study	ScoreCard Example Case	Download Latest
Prostate	79.20y	1.80y	RT00128	Prostate79.20y_2022MAAS_ExampleV2	9/12/2022
Prostate-LN	70.00y (5B)	2.30y	NRG 0103	ProstateLN70.00y_2022MAAS_ExampleV2	12/21/2022
Prostate Bed	48.00y (5B)	20y	NRG 0103 / RADICALS-RT	ProstateBed48.00y_2022MAAS_ExampleV2	9/19/2022
Prostate FLAME	94.50y (5B)	2.70y	FLAME 164	ProstateFLAME94.50y_2022MAAS_ExampleV2	9/19/2022
Prostate SBRT	40.00y (5B)	8.0y	2018 AAMD / RSS Plan Study	ProstateSBRT40.00y_2018AAMD-RSS_PlanStudy	4/15/2016
Breast	50.00y	20y	2014 AAMD / ROR Plan Challenge	WholeBreast50.00y_2014AAMD-ROR_PlanStudy	3/26/2014
Breast	42.560y	2.660y	RT00108	WholeBreast42.560y_2022MAAS_ExampleV2	9/19/2022
Breast-LNs	50.00y	20y	NSABP-RTOG1324/Alliance A222064	ComprehensiveBreast50.00y_2022MAAS_ExampleV2	9/19/2022
Breast-LNs	50.00y	20y	NSABP-RTOG1324	ComprehensiveBreast50.00y_2022MAAS_ExampleV2	9/19/2022
Breast-APBI	30.00y	6.0y	2016 Plan Completion	BreastLateStage30.00y_2016PlanCompletion	2/8/2016
Breast-APBI	30.00y	6.0y	Flowers trial	APBI30.00y_2022MAAS_ExampleV2	9/19/2022
Lung/Mediastinum	50.00y	20y	NRG L034 / Atkins KM 2021	Lung50.00y_2022MAAS_ExampleV2	9/19/2022
Lung	61.20y	1.80y	2021 AAMD Plan Study	Lung61.20y_2021AAMD_PlanStudy	4/26/2021
Lung	69.00y	1.90y	2012 AAMD / ROR Plan Challenge	Lung69.00y_2022MAAS_ExampleV1	5/21/2022
Lung-SBRT	40.00y	20.00y	2014 AAMD / ROR Plan Challenge	Lung40.00y_2014AAMD-ROR_PlanStudy	1/30/2014
Head-and-Neck	70.00y (5B)	20y	NRG 0106	H&N70.00y_2022MAAS_ExampleV2	9/12/2022
Head-and-Neck	70.00y (5B)	20y	2013 Prostate / ROR Plan Challenge	H&N70.00y_2013Prostate-ROR_PlanStudy	10/26/2013
Head-and-Neck	70.00y (5B)	20y	2013 Prostate / ROR Plan Challenge	H&N70.00y_2013Prostate-ROR_PlanStudy	10/26/2013
Esophagus	50.40y (5B)	1.80y	2019 Inter Plan Challenge	Esophagus50.40y_2019InterPlanStudy	12/21/2019
Pancreas	52.00y (5B)	20y	RT00384	Pancreas52.00y_2022MAAS_ExampleV1	5/25/2022
Pancreas	52.00y (5B)	20y	2012 AAMD / ROR Plan Challenge	Pancreas52.00y_2012AAMD_PlanStudy	7/23/2012
Pancreas SBRT	40.00y	8.0y	UK Consensus	PancreasSBRT40.00y_2022MAAS_ExampleV2	9/19/2022
Pancreas-Liver Meta	50.00y (5B)	10.0y	2022 AAMD Plan Study	PancreasLiverMeta50.00y_2022AAMD_PlanStudy	4/25/2022
Liver SBRT	50.00y	10.0y	UK Consensus/RT00112	LiverSBRT50.00y_2022MAAS_ExampleV2	9/19/2022
Spine SBRT	25.00y	25.0y	2011 ROR / Prostate Plan Study	SpineSBRT25.00y_2011ROR_PlanStudy	1/19/2011
Brain	40.00y	20y	RT00825	Brain40.00y_2022MAAS_ExampleV2	9/21/2022
Brain-HSWB	30.00y	30y	NRG-C001	HCSWB30.00y_2022MAAS_ExampleV1	3/28/2022
Anus-LN	50.40y (5B)	1.80y	2013 AAMD / ROR Plan Challenge	Anus50.40y_2013AAMD-ROR_PlanStudy	5/25/2022
Anus	50.40y (5B)	1.80y	RT00280	Anus50.40y_2022MAAS_ExampleV1	6/22/2021
OVN	45Gy	1.80y	NRG 0106	OVN45.00y_2022MAAS_ExampleV2	8/2/2022
TMJ Female	50.00y (5B)	20y	2018 Prostate Plan Challenge	TMJFemale50.00y_2018ProstatePlanStudy	1/24/2018
TMJ Male	20.00y	20y	2018 Prostate Plan Challenge	TMJMale20.00y_2022MAAS_ExampleV1	5/25/2022
Prostate Fossa-LN	48.00y (5B)	20y	2011 AAMD / ROR Plan Challenge	ProstateFossaLN48.00y_2011AAMD-ROR_PlanChallenge	9/19/2011
Liver SBRT	50.00y	10.0y	RT00112	LiverSBRT50.00y_RT00112_ExampleV1	6/9/2022
Breast	26.00y	5.20y	FAST-Forward Trial	WholeBreast26.00y_2022MAAS_ExampleV2	9/19/2022
Breast-LN	57.00y (5B)	2.280y		LiverSBRT57.00y_2022MAAS_ExampleV1	12/09/2022
Breast-LN	57.00y (5B)	2.280y		RtBreast57.00y_5B_2022MAAS_ExampleV1	12/09/2022
Breast-LN	57.00y (5B)	2.280y		Breast57.00y_5B_2022MAAS_ExampleV1	12/09/2022
Brain-HSWB	20.00y	4.0y	CCTO CE.7	HCSWB20.00y_2022MAAS_ExampleV1	12/20/2022

84

Resources

Scorecard Tool

Open-source, ESAPI plan scorecard tool that was used in this project. Create your own dosimetric scorecards and score plans. Instructions for download and usage are available.

<https://github.com/Varian-MedicalAffairsAppliedSolutions/MAAS-PlanScoreCard>

The screenshot displays the GitHub repository page for 'Varian-MedicalAffairsAppliedSolutions / MAAS-PlanScoreCard'. The repository is public and has 253 commits. The 'Releases' section is highlighted with a red box, showing the latest release: 'v16.1-PlanScoreCard-v3.17.12-0...'. Below the release, there is a 'Packages' section and a 'Contributors' section.

85

Acknowledgments

Anthony Magliari
 Lesley Rosa
 Elena Czeizler
 Matt Schmidt
 Jason Paisley
 Vanessa Magliari



Questions?

86