

Hippocampal Sparing Whole Brain: RapidPlan Model Validation Process



Robert Foster, MS DABR CMD; Anthony Magliari, MS CMD; Ryan Clark, MS CMD

Background

Knowledge Based Planning (KBP) is a solution that can expedite the treatment planning process while improving dosimetric standardization across treatment plans. RapidPlan™ is a commercial implementation of KBP (Varian Medical Systems, Palo Alto, CA). In 2016, a Hippocampal Sparing Whole Brain Model Version 1.0 (HSWBv1) was made publicly available.² From this model, we developed a new HSWB KBP model (HSWBv2).¹ This new model placed further emphasis on improving hippocampal sparing, homogeneity, and conformity. A scorecard was used to quantify and compare the relative dosimetric performance (score) for each method between HSWBv1 and HSWBv2.³

Figure 1: Example Version 2.0 scorecard; validation case score 131.4/142




Methods

Upon completion, HSWBv2 was then validated for use with various arc geometries, treatment machines, beam energies, calculation, and optimization methods. HSWBv2 utilized a training set of 42 cases. Five additional validation cases, not included in the original model training set, were used for further evaluating the plan quality between HSWBv1 and HSWBv2. Different treatment delivery systems (Halcyon, Truebeam, and C-series), beam arrangements (coplanar/non-coplanar, HyperArc™, and 3-arc/4-arc VMAT), optimization settings (convergence modes), and calculations algorithms (AcurosXB/AAA) to validate that the HSWBv2.


Table 1: Validation results for different beam arrangements (6X-FFF, AcurosXB v17, extended convergence mode, MR3 return, 2x Intermediate dose)

Patient	Validation Cases HSWB MR3+MR3											
	Halcyon				TrueBeam							
	4 Arcs (Coplanar)		3 Arcs (Coplanar)		4 Arcs (Non-Coplanar)		3 Arcs (Coplanar)		4 Arcs (Coplanar)		HyperArc (Non-Coplanar)	
	V1.0	V2.0	V1.0	V2.0	V1.0	V2.0	V1.0	V2.0	V1.0	V2.0	V1.0	V2.0
36	120.38	132.64	118.46	131.55	121.41	128.72	112.73	126.39	113.35	127.36	118.23	129.4
37	122.87	132.91	119.37	131.92	123.41	129.42	116.91	128.69	117.4	129.46	120	130.1
39	120.57	132.97	119.24	132.89	122.69	129.36	114.17	126.79	113.98	126.3	121.45	130.4
40	119.8	133.56	120.33	132.72	118.08	129.89	112.89	125.9	112.85	126.65	119.54	130.31
41	118.81	132.94	120.41	131.98	114.65	128.01	113.92	126.76	112.95	126.19	122.57	129.44
Average	120.49	133.00	119.56	132.21	120.05	129.08	114.12	126.91	114.11	127.19	120.36	129.93



Different beam arrangements were tested on both Halcyon (SX2MLC) and TrueBeam (M120 MLC). Both the three and four arc Halcyon cases utilized collimator of 315°, (0°), 45°, and 90°. The four arcs non-coplanar Truebeam consisted of two full arcs (0° couch) and two vertex arcs (90° couch) rotating from 180° (PA) -> 5° (from AP) with all arcs utilized a 315°/45° collimator. The 3 arc coplanar plan utilized 315°, 45°, and 90° collimator. The 4 arc plans utilized the same collimator as the 3 arcs except 90° split X jaw superior/inferior to the hippocampus. HyperArc utilized the full 4 arc arrangement. This beam arrangement validation resulted in higher scores across all variations for HSWBv2 against HSWBv1.

Figure 2: Example DVH showing dose difference in hippocampal sparing and PTV coverage between Halcyon (4 arcs), Truebeam (3 arcs), and HyperArc plans.



Results

Convergence mode testing showed that longer optimization time and earlier/multiple returns correlated to higher plan scores. The use of AcurosXB resulted in improved scoring against AAA as well.

Figure 3: Dose calculation and beam energy (Eclipse v17 algorithms, extended convergence mode, MR3 return for intermediate dose)





Figure 4: Convergence Mode: Off, On, Extended (Halcyon 4 arc, 6X-FFF, AcurosXB v17)



Conclusions

The HSWBv2 model validation data shows improved results over HSWBv1. The refined scorecard metrics aided in better quantifying desired improvements according to community feedback from using HSWBv1. Testing HSWBv2 across multiple delivery systems, beam geometries, beam energies, optimization processing, and algorithms provides additional confidence in model quality and future use expectations. This KBP model is available for public use along with further information on how to implement the model in a clinical description document.⁴

Reference

- (1) Liu, Hefei et al. "RapidPlan hippocampal sparing whole brain model version 2-how far can we reduce the dose?" *Medical dosimetry: official journal of the American Association of Medical Dosimetrists*, S0958-3947(22)00039-5. 2 May. 2022, doi:10.1016/j.meddos.2022.04.003.
- (2) Magliari, Vanessa et al. "Hippocampal sparing whole brain: Rapid Plan Model Following the NRG-CC001 Protocol". AAMD conf poster present. 2017; (Published online available at:) <https://medicallaffairs.varian.com/download/PosterPresentationAAMD2017RapidPlanHCSWB.pdf>
- (3) <https://github.com/Varian-innovation-center>
- (4) <https://medicallaffairs.varian.com/wholebrain-hippocampalsparing-vmat2>